

Self-dissolving paper technology is ideal for fast-paced Operating Rooms



Drop in before or after fluid is collected

Drop-In Pacs

NEW! Advanced Self-Dissolving Paper Technology

The **Drop-In Pacs** from Safetec are a newly designed medical solidifier with an advanced self-dissolving technology to instantly absorb and solidify fluid waste. Containing the same effective super absorbent polymers as used in Red Z and Green-Z bottles, these slender pouches are now made with a faster acting, innovative water-soluble paper.

For use in: kick buckets, operating rooms, red biohazard bags, suction canisters and wound care therapy



* Red Z deodorizes spills (contains chlorine).

Improve OR Turnover and Worker Safety

- Specifically designed to fit in porthole of canister
- Easy and ready to use premeasured pouches
- Faster acting with new water-soluble paper technology
- Drop-In Pacs are eco friendly and reduce extra disposal costs that bottles require
- Contents are non-toxic



Directions for Use

For Suction Canister: Insert entire drop-in pac through port opening of suction canister (before or after surgical procedure). Allow sufficient time for encapsulation before handling.

For Biohazard Bag or Kick Bucket: Place drop-in pac in the bottom of a biohazard bag or other container to control fluid before or after spills occur. Clean and disinfect contaminated areas after use. Dispose of infectious waste in accordance with all local, state and federal regulations.

Drop-In Pac Absorption Rates

Item#	Description	Time to Open	Time to Absorb	Capacity
41108	Red Z® Drop-In (20 gr)	10–25 sec	1–3 min	750–900cc
41175	Red Z® Drop-In (40 gr)	10–25 sec	2–6 min	1500–1800cc
42008	Green-Z® Drop-In (20 gr)	10–25 sec	1–3 min	750–900cc
42175	Green-Z® Drop-In (40 gr)	10–25 sec	2–6 min	1500–1800cc

* all packets tested in 0.9% saline solution in a temperature range of 32–37°C



OSHA Compliance 1910.1030 – Bloodborne Pathogens Standard

Universal Precautions is an approach to infection control. According to the concept of Universal Precautions, all human blood and certain human body fluids are treated as if known to be infectious.