



Examine Real Quality™

BLOOD DRAWING CHAIRS



SPECIALTY CHAIRS

ABOUT CLINTON

No other manufacturer of medical equipment, knows blood chairs exactly like the medical equipment experts at Clinton Industries. In fact, we take pride in being among the industry's leaders in the manufacture and sale of blood drawing chairs. Maybe it's our broad line of products. Maybe it's our two decade reputation for quality and integrity. Or perhaps, it's our knowledgeable and friendly customer service that has made Clinton a prime industry source for blood drawing chairs. Discover what much of the industry already knows and *examine real quality*—Clinton Industries, for blood drawing chairs and all of your medical equipment needs.

See the entire line of Clinton product catalogs

- *Medical Tables & Cabinets*
- *Pediatric Medical Equipment*
- *Medical Seating*
- *Medical Accessories*
- *Physical Therapy Equipment*
- *Quick Ship Tables*



6000X

Extra-Wide Blood Drawing Chair

- Accommodates larger patient or child with parent
- Padded & upholstered seat and arms
- Adjustable height arms 29"–34"
- Two pivot arms
- All welded heavy-duty frame

Seat—34" wide x 18" deep x 20" high
Overall—46" wide x 26" deep



6099X

Extra-Wide and Tall Blood Drawing Chair

- Same as above with extra seat height—26" and footrest
- Reduces constant bending for staff
- Extra height makes sitting and rising easier for older patients



6077

Infant Blood Drawing Station

- All welded stainless steel frame
- Top with two rims
- Two safety straps
- Plastic work tray
- Removable middle padded section

Top—20" wide x 32" long


Overall—20 1/2" wide x 39" long x 39 1/2" high



6088

Wheel Chair Blood Drawing Station

- Adjustable height 26"–39"
- Adjustable depth 24"–36"
- Padded armrest 6" x 24"
- Mounts to wall and floor
- Plastic work tray (optional)
- Arm folds up for easy storage when not in use



DURA LINE

Our best line of blood drawing chairs, Dura Line chairs are extra strong and extra comfortable.

FEATURES

- All heavy-gauge steel, welded frame
- Powder coated frame
- Comfortable 2" padded seat and backrest
- 360° swing arm for easy alignment
- Fully adjustable height and width arm
- Optional padded flip arm available
- 20 colors of vinyl



6020

Blood Drawing Chair with Drawer

- Padded & upholstered seat and arms
- One piece plastic drawer
- Adjustable height arms 27"–33"
- Pivot armrest

Seat—17" wide x 16" deep x 20" high
Overall—43" wide x 26" deep



6010

Blood Drawing Chair

- Padded & upholstered seat and arms
- Adjustable height arms 27"–33"
- Pivot armrest

Seat—17" wide x 16" deep x 20" high

Overall—34" wide x 26" deep



6022

Extra-Tall Blood Drawing Chair with Drawer

- Reduces constant bending for staff
- Extra height makes sitting and rising easier for older patients
- Padded & upholstered seat and arms
- Adjustable height arms 34"–40"
- Easy clean one piece plastic drawer
- Patient foot rest
- Pivot armrest

Seat—17" wide x 16" deep x 26" high

Overall—43" wide x 26" deep



6011

Extra-Tall Blood Drawing Chair

- Same as above without drawer

Overall—34" wide x 26" deep

OPTIONS



680 Optional Adjustable Length Flip-Arm

- Adjustable arm height
- Adjustable depth for larger patients
- Fully padded
- Fits all model chairs

6" wide x 23" long

STYLE LINE

The Style Line offers just that—a stylish looking alternative to traditional designs.

FEATURES

- All heavy-gauge steel, welded frame
- Powder coated frame
- Comfortable 2" Padded seat
- 360° swing arm for easy alignment
- Fully adjustable height and width arm
- Optional padded flip arm available
- 20 colors of vinyl



64929

Blood Chair with Drawer

- Padded & upholstered seat and arm pads
- Adjustable height arms 28"–33"
- Pivot armrest
- One piece plastic drawer
- Adjustable arm width 19"–24"

Seat—17" wide x 16" deep x 20" high
Overall—42" wide x 26" deep



64929-B

Blood Chair with Drawer

- Padded & upholstered seat, armrests and back
- Adjustable height arms 28"–33"
- Pivot armrest
- White chair frame
- One piece plastic drawer
- Adjustable arm width 19"–24"

Seat—17" wide x 16" deep x 20" high
Overall—42" wide x 26" deep



64950

Blood Drawing Chair

- Padded & upholstered seat and armrests
- Adjustable height arms 28"–33"
- Pivot armrest
- White chair frame
- Adjustable arm width 19"–24"

Seat—17" wide x 16" deep x 20" high
Overall—32" wide x 26" deep



64950-B

Blood Drawing Chair

- Padded & upholstered seat, back, and armrests
- Adjustable height arms 28"–33"
- Pivot armrest
- White chair frame
- Adjustable arm width 19"–24"

Seat—17" wide x 16" deep x 20" high
Overall—32" wide x 26" deep

VALUE SERIES

Economical and quality describe the value series of chairs. Check these Clinton features and you'll see why.

FEATURES


- All steel, welded frame
- Fully adjustable height and arm position
- Optional padded flip arm available
- 360° swing arm for easy alignment



6060-P
Blood Drawing Chair
• One piece plastic seat
• Adjustable height arms 27"–33"
• Laminate pivot arm rest
Seat — 17" wide x 16" deep x 20" high
Overall — 34" wide x 26" deep



6029-P
Blood Drawing Chair
• One piece plastic seat
• Adjustable height arms 27"–33"
• Laminate pivot armrest
• Easy-clean one-piece drawer
Seat — 17" wide x 16" deep x 20" high
Overall — 43" wide x 26" deep




6011-P
Extra-Tall Blood Drawing Chair
• Reduces constant bending
• One piece plastic seat
• Adjustable height arms 34"–40"
• Laminate pivot armrest
• Patient foot rest
Seat—17" wide x 16" deep x 26" high
Overall—34" wide x 26" deep



6022-P
Extra-Tall Blood Drawing Chair
• Reduces constant bending
• One piece plastic seat
• Adjustable height arms 34"–40"
• Laminate pivot armrest
• Patient foot rest
• Easy-clean one-piece drawer
Seat—17" wide x 16" deep x 26" high
Overall—43" wide x 26" deep

OPTIONS



671
Optional Laminate Flip-Arm
• Fits on all model chairs
• Arm is adjustable in height and depth
• Arm does not pivot. It remains horizontal in front of patient
• Easily flips out of way
• Full-length armrest allows for proper positioning of either arm
6" wide x 23" long



674
Optional Padded Armrest
• Same as standard laminate armrest with the added comfort of upholstery.