

SEATING COMPARISON CHART

Which Clinton stool is right for you? Compare the key features of all of our fine seating choices quickly and easily with the chart below. Then place your order for the stools that fit your needs.

	Model No.	Casters	Height Range	Base	Adjustment	Rests
SPECIALTY	2155-SP	Nylon	18"-23 ¹ / ₂ "	Aluminum	Pneumatic	Crescent Armrest
	2145-SP	Nylon	19 ¹ / ₂ "-24 ¹ / ₂ "	Aluminum	Pneumatic	Crescent Armrest
	2145-W	Nylon	18"-23"	Aluminum	Pneumatic	Backrest
	2145	Nylon	19 ¹ / ₂ "-24 ¹ / ₂ "	Aluminum	Pneumatic	None
EPIC	2188-W	Nylon	21" - 29"	Aluminum	Pneumatic	Backrest & Armrests
	2175-W	Nylon	16 ¹ / ₂ "-22"	Aluminum	Pneumatic	Backrest
	2167-W	Nylon	17"-22 ¹ / ₂ "	Aluminum	Pneumatic	Backrest
	2166-W	Nylon	21"-29"	Aluminum	Pneumatic	Backrest
	2165-W	Nylon	17"-22 ¹ / ₂ "	Aluminum	Pneumatic	Backrest & Armrests
SELECT	2166	Nylon	22 ¹ / ₂ "-30"	Aluminum	Pneumatic	Backrest
	2165	Nylon	22 ¹ / ₂ "-30"	Aluminum	Pneumatic	Backrest
	2155	Nylon	18"-23 ¹ / ₂ "	Aluminum	Pneumatic	Optional Backrest
	2150	Nylon	17 ¹ / ₄ "-21 ¹ / ₂ "	Aluminum	Thread	Optional Backrest
STYLE LINE	2135-S	Nylon	18"-23"	Slate Nylon	Pneumatic	Optional Backrest
	2135-P	Nylon	18"-23"	Putty Nylon	Pneumatic	Optional Backrest
	2135	Nylon	18"-23"	Black Nylon	Pneumatic	Optional Backrest
	2130	Nylon	18"-22 ¹ / ₄ "	Black Nylon	Thread	Optional Backrest
CLASSICS	2102-CC	Nylon	18 ¹ / ₂ "-25 ¹ / ₂ "	Chrome	Thread	None
	2100-CC	Nylon	18 ¹ / ₂ "-25 ¹ / ₂ "	Chrome	Thread	None
	2102	Rubber	18 ¹ / ₂ "-25 ¹ / ₂ "	Chrome	Thread	None
	2100	Rubber	18 ¹ / ₂ "-25 ¹ / ₂ "	Chrome	Thread	None
	2192	Rubber	23"-29 ¹ / ₂ "	Chrome	Thread	None