

DURA LINE

Our best line of blood drawing chairs, Dura Line chairs are extra strong and extra comfortable.

FEATURES

- All heavy-gauge steel, welded frame
- Powder coated frame
- Comfortable 2" padded seat and backrest
- 360° swing arm for easy alignment
- Fully adjustable height and width arm
- Optional padded flip arm available
- 20 colors of vinyl



6020

Blood Drawing Chair with Drawer

- Padded & upholstered seat and arms
- One piece plastic drawer
- Adjustable height arms 27"–33"
- Pivot armrest

Seat—17" wide x 16" deep x 20" high
Overall—43" wide x 26" deep



6010

Blood Drawing Chair

- Padded & upholstered seat and arms
- Adjustable height arms 27"–33"
- Pivot armrest

Seat—17" wide x 16" deep x 20" high

Overall—34" wide x 26" deep



6022

Extra-Tall Blood Drawing Chair with Drawer

- Reduces constant bending for staff
- Extra height makes sitting and rising easier for older patients
- Padded & upholstered seat and arms
- Adjustable height arms 34"–40"
- Easy clean one piece plastic drawer
- Patient foot rest
- Pivot armrest

Seat—17" wide x 16" deep x 26" high
Overall—43" wide x 26" deep




6011

Extra-Tall Blood Drawing Chair

- Same as above without drawer

Overall—34" wide x 26" deep

OPTIONS



680 Optional Adjustable Length Flip-Arm

- Adjustable arm height
- Adjustable depth for larger patients
- Fully padded
- Fits all model chairs

6" wide x 23" long