


EXERCISE PRODUCTS



L-2630

L-2630

Ramp and Curb Training Set

An adjustable training unit in which curbs can be removed or sequenced according to training requirements. Curbs are stored beneath platform.

Curbs: 32" x 32" x 2" high, 36" x 36" x 4" high,

40" x 40" x 6" high

Ramp: 84 x 44 x 8 (L x W x H).

Incline extends from 40° to 154.95584°



LBP-600



ASW-413



AWP-420



AT-440



A-4020

AT-440

Over Door Traction Unit

- Suitable for home or office traction
- Support elements, spreader bar, nylon cord and graduated water bag

A-4020

Shoulder Ladder

- 36 steps for upper extremity exercise
- Wall mounted
- Wood with lacquer finish



Examine Real Quality™

545 S. Pine St., York, PA 17403

800-441-9131 • 717-848-3519

Fax 717-843-5871

www.clinton-ind.com