

# FREE STANDING CABINETS



	Length	Height	Depth
<b>8821</b>	24"	35"	18"

### Treatment Cabinet with 2 Doors & 1 Drawer

- 1 adjustable shelf
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction



	Length	Height	Depth
<b>8822</b>	24"	35"	18"

### Treatment Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction



	Length	Height	Depth
<b>8805</b>	24"	35"	18"

### Treatment Cabinet with 5 Drawers

- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction



	Length	Height	Depth
<b>8922</b>	24"	34"	18"

### Mobile Treatment Cabinet with 2 Doors and 2 Drawers

- Mounted on dual wheel swivel casters
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction



	Length	Height	Depth
<b>8920</b>	22"	33"	17"

### Mobile Equipment Cart with 2 Drawers

- Mounted on a stable 5-leg black nylon base
- 2 Euro-style drawers
- Top rimmed on 3 sides
- All laminate cabinet



	Length	Height	Depth
<b>8940</b>	22"	33"	17"

### Mobile Equipment Cart with 1 Drawer and 1 Shelf

- Mounted on a stable 5-leg black nylon base
- 1 Euro-style drawer
- Top and back have wire feed cut-outs
- Top rimmed on 3 sides
- All laminate cabinet



	Length	Height	Depth
<b>089</b>	17½"	5"	6¾"

### Optional Basket

- Wire basket with wire epoxy finish