

MAT PLATFORMS

STANDARD PLATFORM FEATURES

- Six solid oak legs held in place with three bolts each
- Extra center support
- Ultra firm 2" thick padding
- Custom bumper strips protects vinyl top
- Frame secured by steel angle
- 18" standard height with 21" height option
- 20 colors of upholstery



MO-240
Upholstered Mat Platform

- MO-240-46** 4' x 6'
- MO-240-47** 4' x 7'
- MO-240-57** 5' x 7'
- MO-240-68** 6' x 8'



CH-240
Electric Lift Mat Platform
Power mat platform with upholstered top

- CH-240-46** 4' x 6'
- CH-240-47** 4' x 7'
- CH-240-57** 5' x 7'
- CH-240-68** 6' x 8'




MO-242
Folding Mat Platform
Projects only to 8" when folded against the wall

- MO-2430-47** 4' x 7'
- MO-2430-57** 5' x 7'

POWER PLATFORM FEATURES

- All welded steel frame
- 700 lb. lifting capacity
- Foot pedal control frees practitioner's hands
- Motor enclosed in acrylic panels for quiet operation
- Single motor drive
- Actuator lifts from low to high in 15 seconds
- Overload protection
- 110 volt, 60 Hz, 2.7 amps
- Bumper guards protect vinyl
- 20 colors of upholstery



CH-240
Power Mat Platform
Mat attaches with Velcro® strips

- CH-242-47** 4' x 7'
- CH-242-57** 5' x 7'
- CH-242-68** 6' x 8'



Due to the limitations of the printing process, actual colors may vary.