

PARALLEL BARS



FEATURES

- One piece stainless steel handrails with bumper ends
- Uprights are ultra-heavy telescoping steel tube
- Adjustable height and width
- Height controlled by spring loaded plungers for easy adjustment
- Any play is eliminated by a second knob at each upright
- Platforms have durable satin topcoat and solid hardwood ends
- Safety treads at both ends on platform models
- Options available to customizable your equipment



PFL-3410

Floor Mount Parallel Bars

Floor mounted, 10' long stainless steel handrails are adjustable in both height and width.

Height 26" to 39"

Width adjusts out 6" on each side.

PP-2012

Same as above but 12' long.

PP-2007


Platform Mounted Parallel Bars

7' long stainless steel handrails

Height: 26" to 39" • Width: 19" to 28"

Rail length: 7'

Note: For child rails & abduction board options please call.



PP-2010

Same as above but 10' long

PP-2012

Same as above but 12' long.

PF-33W

Wall Folding Parallel Bars

The best choice for parallel bars in a limited space situation. Seven feet long, rounded handrails fold up against the wall when not in use. Only requires 14" of floor space when not in use.

Height adjusts from 28" to 42" • Length is 7'



PF-30A

Folding Parallel Bars

7' long stainless steel handrails with easy storage option when not in use. Unique stabilizer bars stiffen the rails. Heavy gauge steel uprights with spring plungers and sure lock, positioning knobs allow for sure, easy height change.

Height: 22' to 35" • Width: 16" to 24" • Length: 7'

PF-30X

Same as above but 10' long.