



Examine Real Quality™

PHYSICAL THERAPY EQUIPMENT



GLASS MIRRORS

ABOUT CLINTON

Getting back on the road to recovery is hard work for the patient and the therapist. That's why Clinton makes dependable, quality constructed physical therapy equipment. We want our equipment to work just as hard as the patients and medical professionals who use them day in and day out. Discover what much of the medical profession already knows—*examine real quality*—Clinton Industries, for physical therapy equipment and all of your medical equipment needs.

See the entire line of Clinton product catalogs

- Medical Tables & Cabinets
- Pediatric Medical Equipment
- Blood Drawing Chairs
- Medical Seating
- Medical Accessories
- Quick Ship Tables

These fine Signature Series mirrors feature ANSI safety backed mirror glass framed in the classic beauty of natural solid oak—crafted the Clinton way for long lasting service.



2201-9MO
Mobile Mirror
27" wide x 72" high



2221-8MO



2200-MO



2230-MO

2221-8MO
Child's Mobile Mirror
19" wide x 58" high

2200-MO
Wall Mount Mirror
27" wide x 69" high

2230-MO
Tri-Panel Mirror
60" wide x 72" high

PARALLEL BARS



FEATURES

- One piece stainless steel handrails with bumper ends
- Uprights are ultra-heavy telescoping steel tube
- Adjustable height and width
- Height controlled by spring loaded plungers for easy adjustment
- Any play is eliminated by a second knob at each upright
- Platforms have durable satin topcoat and solid hardwood ends
- Safety treads at both ends on platform models
- Options available to customizable your equipment



PFL-3410

Floor Mount Parallel Bars

Floor mounted, 10' long stainless steel handrails are adjustable in both height and width.

Height 26" to 39"

Width adjusts out 6" on each side.

PFL-3412

Same as above but 12' long.

PP-2007


Platform Mounted Parallel Bars

7' long stainless steel handrails

Height: 26" to 39" • Width: 19" to 28"

Rail length: 7'

Note: For child rails & abduction board options please call.



PP-2010

Same as above but 10' long

PP-2012

Same as above but 12' long.

PF-33W

Wall Folding Parallel Bars

The best choice for parallel bars in a limited space situation.

Seven feet long, rounded handrails fold up against the wall when not in use. Only requires 14" of floor space when not in use.

Height adjusts from 28" to 42" • Length is 7'



PF-30A

Folding Parallel Bars

7' long stainless steel handrails with easy storage option when not in use. Unique stabilizer bars stiffen the rails. Heavy gauge steel uprights with spring plungers and sure lock, positioning knobs allow for sure, easy height change.

Height: 22' to 35" • Width: 16" to 24" • Length: 7'

PF-30X

Same as above but 10' long.

STAIRCASES



S-5080-30



Above: Straight Position

At Left: Corner Position



S-5000-30



S-5010-30



A-500-30

TRAINING STAIRS *with steel uprights and rounded wood handrails*

S-5080-30

Staircase converts from straight to corner arrangement

Platform 30" long to 30" wide

Four Steps: 30" wide x 6"high/ Six Steps: 30" wide x 4"high

S-5000-30

Small Staircase

Small straight staircase for limited space

Platform 24" long to 30" wide

Two Steps: 30" wide x 6"high/ Three Steps: 30" wide x 4"high

Note: Shipped assembled

S-5010-30

One-Side Staircase

Staircase with platform and three enclosed sides

Platform 24" long to 30" wide • Four Steps: 30" wide x 6"high

Note: Shipped assembled

TRAINING STAIRS *with adjustable height stainless steel handrails*

A-501-30

Adjustable One-Side Staircase

Adjusts to all ages from pediatric to adult in 1 1/2" increments

Platform 24" long to 30" wide

Four Steps: 30" wide x 6"high

A-500-30

Adjustable Small Staircase

Similar to A-501-30, but patient can ascend and descend without turning

Platform 24" long to 30" wide

Two Steps: 30" wide x 6"high

Three Steps: 30" wide x 4"high



A-501-30

MAT PLATFORMS

STANDARD PLATFORM FEATURES

- Six solid oak legs held in place with three bolts each
- Extra center support
- Ultra firm 2" thick padding
- Custom bumper strips protects vinyl top
- Frame secured by steel angle
- 18" standard height with 21" height option
- 20 colors of upholstery



MO-240
Upholstered Mat Platform

- MO-240-46** 4' x 6'
- MO-240-47** 4' x 7'
- MO-240-57** 5' x 7'
- MO-240-68** 6' x 8'



CH-240
Electric Lift Mat Platform
Power mat platform with upholstered top

- CH-240-46** 4' x 6'
- CH-240-47** 4' x 7'
- CH-240-57** 5' x 7'
- CH-240-68** 6' x 8'




MO-242
Folding Mat Platform
Projects only to 8" when folded against the wall

- MO-2430-47** 4' x 7'
- MO-2430-57** 5' x 7'

POWER PLATFORM FEATURES

- All welded steel frame
- 700 lb. lifting capacity
- Foot pedal control frees practitioner's hands
- Motor enclosed in acrylic panels for quiet operation
- Single motor drive
- Actuator lifts from low to high in 15 seconds
- Overload protection
- 110 volt, 60 Hz, 2.7 amps
- Bumper guards protect vinyl
- 20 colors of upholstery



CH-240
Power Mat Platform
Mat attaches with Velcro® strips

- CH-242-47** 4' x 7'
- CH-242-57** 5' x 7'
- CH-242-68** 6' x 8'



Due to the limitations of the printing process, actual colors may vary.

EXERCISE PRODUCTS

L-2630

Ramp and Curb Training Set

An adjustable training unit in which curbs can be removed or sequenced according to training requirements. Curbs are stored beneath platform.

Curbs: 32" x 32" x 2" high, 36" x 36" x 4" high,

40" x 40" x 6" high

Ramp: 84 x 44 x 8 (L x W x H).

Incline extends from ³