

STAINLESS STEEL PRODUCTS



| | Length | Height | Width |
|---------------|---------|--------|---------|
| SS-150 | 14 1/4" | 9" | 11 1/4" |

Stainless Steel Step Stool with Hand Rail

- All welded construction
- Rubber tread top for safety
- 7/8" diameter handrail
- Rubber feet

| | Length | Height | Width |
|---------------|---------|--------|---------|
| SS-140 | 14 1/4" | 9" | 11 1/4" |

Stainless Steel Step Stool

- All welded construction
- Rubber tread top for safety
- Rubber feet



Diameter Total Height

| | | |
|---------------|---------|-----|
| SS-160 | 14 1/2" | 12" |
|---------------|---------|-----|

Stainless Steel Kick Bucket and Frame

- Heavy duty stainless steel frame
- 360° bumper on frame
- 2" rubber wheel ball bearing casters
- 13 quart stainless steel bucket

Diameter

| | |
|---------------|---------|
| SS-166 | 14 1/2" |
|---------------|---------|

Stainless Steel Kick Bucket Frame Only



Seat Diameter Height Range

| | | |
|----------------|-----|-------------|
| SS-2142 | 14" | 25"–29 1/2" |
|----------------|-----|-------------|

Stainless Steel Stool

- All welded frame
- Smooth machine screw height adjustment
- 2" rubber wheel ball bearing casters
- 4" Thick seat for comfort
- Round foot ring
- Easy clean upholstery

Seat Diameter Height Range

| | | |
|-------------------|-----|-------------|
| SS-2142-21 | 14" | 25"–29 1/2" |
|-------------------|-----|-------------|

Stainless Steel Stool with Backrest

- Same as SS-2142 with padded backrest

Base Diameter Height Range

| | | |
|--------------|-----|----------|
| SS-31 | 22" | 60"–103" |
|--------------|-----|----------|

Stainless Steel IV Pole

- Stainless steel construction
- 2 ram's horn hooks
- 2" rubber wheel ball bearing casters

Base Diameter Height Range

| | | |
|--------------|-----------|---------|
| SS-21 | 18" x 22" | 31"–50" |
|--------------|-----------|---------|

Stainless Steel Instrument Stand

- Stainless steel frame
- Removable stainless steel tray
- U-shaped base
- 2" rubber wheel ball bearing casters
- Single knob adjustment

NOTE: The above stainless steel items are not to be used near an M.R.I.