

# STYLE LINE

The Style Line offers just that—a stylish looking alternative to traditional designs.

## FEATURES

- All heavy-gauge steel, welded frame
- Powder coated frame
- Comfortable 2" Padded seat
- 360° swing arm for easy alignment
- Fully adjustable height and width arm
- Optional padded flip arm available
- 20 colors of vinyl



### 64929

#### Blood Chair with Drawer

- Padded & upholstered seat and arm pads
- Adjustable height arms 28"–33"
- Pivot armrest
- One piece plastic drawer
- Adjustable arm width 19"–24"

**Seat—17" wide x 16" deep x 20" high**  
**Overall—42" wide x 26" deep**



### 64929-B

#### Blood Chair with Drawer

- Padded & upholstered seat, armrests and back
- Adjustable height arms 28"–33"
- Pivot armrest
- White chair frame
- One piece plastic drawer
- Adjustable arm width 19"–24"

**Seat—17" wide x 16" deep x 20" high**  
**Overall—42" wide x 26" deep**



### 64950

#### Blood Drawing Chair

- Padded & upholstered seat and armrests
- Adjustable height arms 28"–33"
- Pivot armrest
- White chair frame
- Adjustable arm width 19"–24"

**Seat—17" wide x 16" deep x 20" high**  
**Overall—32" wide x 26" deep**



### 64950-B

#### Blood Drawing Chair

- Padded & upholstered seat, back, and armrests
- Adjustable height arms 28"–33"
- Pivot armrest
- White chair frame
- Adjustable arm width 19"–24"

**Seat—17" wide x 16" deep x 20" high**  
**Overall—32" wide x 26" deep**