

# SUPPORT ACCESSORIES



	Length	Width
<b>TS-170 Top</b>	15"	31 1/2"
<b>Base</b>	16"	26 1/2"
<b>Height Range</b>	28"—43"	

### Over Bed Table

- Walnut laminated top
- Chrome plated base
- Dual wheel nylon casters
- Single hand spring loaded height adjustment
- Shipped KD
- Weight capacity 50 lbs.

	Length	Width
<b>TS-175 Top</b>	15"	31 1/2"
<b>Base</b>	16"	26 1/2"
<b>Height Range</b>	28"—43"	

### Over Bed Table

- Gray laminated top
- Chrome plated base
- Dual wheel nylon casters
- Single hand spring loaded height adjustment
- Shipped KD
- Weight capacity 50 lbs.



	Length	Height	Depth
<b>T-32</b>	11 1/2"	18 1/4"	12 1/4"

### 32 Quart Waste Can

- White steel shell
- Foot pedal operated lid
- Removable plastic inner bucket

	Length	Height	Width
<b>T-6142</b>	14 1/4"	9"	11 1/4"

### Bariatrics Foot Stool

- Chrome-plated base
- Reinforced "X" base construction
- Rubber tread top for safety
- Rubber feet
- 600 lbs. capacity



	Height Range
<b>T-10</b>	50"—74"

### Gooseneck Lamp

- UL & CSA approved
- Grounded three-prong, 12' cord
- Chrome finish
- Inner/outer "cool shade" system
- Weighted chrome base
- Flexible gooseneck arm



	Length	Height	Width
<b>T-50</b>	14 1/4"	9"	11 1/4"

### Step Stool with Hand Rail

- Chrome plated steel
- Rubber tread top for safety
- 7/8" diameter handrail
- Rubber feet



	Length	Height	Width
<b>T-40</b>	14 1/4"	9"	11 1/4"

### Step Stool (one per box)

- Chrome plated steel
- Rubber tread top for safety
- Rubber feet

	Length	Height	Width
<b>T-42</b>	14 1/4"	9"	11 1/4"

### Step Stool (two per box)

- Same as T-40