

VALUE SERIES

Economical and quality describe the value series of chairs. Check these Clinton features and you'll see why.

FEATURES


- All steel, welded frame
- Fully adjustable height and arm position
- Optional padded flip arm available
- 360° swing arm for easy alignment



6060-P
Blood Drawing Chair
• One piece plastic seat
• Adjustable height arms 27"–33"
• Laminate pivot arm rest
Seat — 17" wide x 16" deep x 20" high
Overall — 34" wide x 26" deep



6029-P
Blood Drawing Chair
• One piece plastic seat
• Adjustable height arms 27"–33"
• Laminate pivot armrest
• Easy-clean one-piece drawer
Seat — 17" wide x 16" deep x 20" high
Overall — 43" wide x 26" deep



6011-P
Extra-Tall Blood Drawing Chair
• Reduces constant bending
• One piece plastic seat
• Adjustable height arms 34"–40"
• Laminate pivot armrest
• Patient foot rest
Seat—17" wide x 16" deep x 26" high
Overall—34" wide x 26" deep



6022-P
Extra-Tall Blood Drawing Chair
• Reduces constant bending
• One piece plastic seat
• Adjustable height arms 34"–40"
• Laminate pivot armrest
• Patient foot rest
• Easy-clean one-piece drawer
Seat—17" wide x 16" deep x 26" high
Overall—43" wide x 26" deep

OPTIONS



671
Optional Laminate Flip-Arm
• Fits on all model chairs
• Arm is adjustable in height and depth
• Arm does not pivot. It remains horizontal in front of patient
• Easily flips out of way
• Full-length armrest allows for proper positioning of either arm
6" wide x 23" long



674
Optional Padded Armrest
• Same as standard laminate armrest with the added comfort of upholstery.