



D900

Mini Dopplex® Non-directional Doppler

Product #
D900-P-USA

Description:

The Mini Dopplex® pocket audio Doppler with standard features. When combined with the EZ8 probe, it is ideal for vascular assessment and ABI measurements.

Available Probes: 2, 3, 4, 5, 8, EZ8 and 10 MHz



List Price 2003 \$669.00

Additional Copy:

The New High Sensitivity Mini Dopplex® non-directional Doppler is a value priced pocket unit providing all the standard features of a traditional Doppler. As one of the family of advanced Dopplex units, it benefits from the same leading technology and manufacturing techniques used in more sophisticated models to provide the same level of quality, reliability and performance. It is an ideal model for ABI assessment. Interchangeable probes allow flexibility for vascular and obstetric testing. Each model includes a 9-volt battery, carry case, headphones and gel.

Huntleigh Healthcare Inc.
A Huntleigh Technology Inc. Company
40 Christopher Way, Eatontown, NJ 07724-3327
Tel: (800) 223-1218 (732) 578-9898 Fax: (732) 578-9889