

DFK2 with D900 **Diabetic Foot Kit with D900**

Product#
DFK2-USA-D900

Description:

The Diabetic Foot Kit is a complete package for Peripheral Vascular testing in the diabetic patient. This kit comes with the Mini Dopplex[®] Audio Doppler. The kit includes the following: sphyg, toe and arm cuffs, 10 ga. Monofilament Neuropen, training video, gel and VP8 probe.

Applications include:
Ankle Brachial Index (ABI)
Toe Brachial Index (TBI)
Neuropathic detection

