


Dopplex® Reporter Software Package 3.5

Product # DR3

Dopplex Reporter software package provides a solution to today's patient documentation needs. Complete vascular studies can be undertaken and saved in a patient database. Full page documented printouts can be obtained or results can be networked to other locations. The software is easily installed and does not require modification to your computer.

Requires Windows™ 95, 98, NT4, ME, 2000 or XP



Requirements:

Multi Dopplex® II (MD2) handheld Doppler

ACC76 Serial cable (included with the software)

9 Pin serial male type connector OR USB on PC or Desktop Computer