



FD1+
Fetal Dopplex® Handheld Doppler

Product #
FD1+-P-USA

Description:

Fetal Heart Rate Display
Fixed High Sensitivity 2MHz probe
Advanced auto-correlation fetal heart rate and digital display
LCD Display of fetal heart rate
Waterproof probe as standard



List Price 2003 \$776.00

Additional Copy:

The New High Sensitivity Fetal Dopplex® rate display Doppler offers all the traditional features of an obstetric Doppler combined with an advanced heart rate display capability. The rugged, reliable and waterproof 2MHz fixed probe fetal Dopplex Doppler is ideal for use in all healthcare facilities.

Huntleigh Healthcare L.L.C.
A Huntleigh Technology Inc. Company
40 Christopher Way, Eatontown, NJ 07724-3327
Tel: (800) 223-1218 (732) 578-9898 Fax: (732) 578-9889