

FD2
Fetal Dopplex® II Handheld Doppler

Product #
FD2-P-USA

Description:

- Top of the line Obstetric Doppler
- Interchangeable probes
- New High Sensitivity 2 & 3MHz obstetric probe options
- Unique three fetal heart rate processing modes for optimum performance
- Real time heart rate data output for connection to Dopplex® Printa II
- Event marker- record fetal movement on Dopplex® Printa II



The New High Sensitivity Fetal Dopplex® II bi-directional Doppler is the most advanced handheld unit in the world. It provides 3 levels of signal processing for fetal heart rate along with printout capability using the Dopplex® Printa II. It is available with a choice of 2 or 3 MHz obstetric probes. It can also be used as a vascular Doppler when used with a vascular probe. Each model includes a 9 volt battery, carry case, and gel.