

ProPlus® Digital Handrail Scale

When patients are less mobile our ProPlus™ Stand-On Scale with built-in wheels moves easily to them. Extremely durable and safe, the easy grip handrail and non-skid platform provide added support to unsteady patients.

Features include:

- 180° swivel head for side read
- Keypad entry
- Customizable readout
- Oversize platform
- Computer interface with 270 patient memory

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



Health o meter®
Professional

ProPlus® Clinical Stand-On Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



2101KL

- 180° swivel head for side read
- Motion-sensing weighing technology
- Customizable 1" backlit LCD display
- Displays fractions, lb/oz, decimals in standard; decimals in metric
- User-friendly menu screen
- Clock with date and time
- Keypad Entry, BMI, Manual Tare, Zero, Auto-Off, Auto-Hold, Reweigh, Hold/Release
- Industrial-strength wrap-around handrails
- USB and Serial Ports (2) for computer interface with 270 patient memory
- Updatable firmware
- 6 D-Cell Batteries (6 month battery life) & 120V AC adapter
- Heavy-duty wheels
- 3 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
2101KL	800 lb 360 kg	4 oz/0.2 lb/1/4 lb 0.1 kg	66 lb	26" x 24 1/4"	30" x 28 1/2" x 47 1/4"

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
2101KL	027045642453	1	UPC	50" x 29 3/8" x 9 7/8"	8.34	75 lb

Health o meter®
Professional