

Bariatric Digital Stand-On Scale

This high-capacity stand-on scale provides extra sensitivity and support to patients who are unsteady on their feet. Advanced technology compensates for patient movement.

Features include:

- Motion sensing weighing technology
- Heavy-duty wheels for mobility
- Large 1" LCD display

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



1101KL

Health o meter®
Professional

Electronic Clinical Stand-On Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



1101KL

- Easy-grip handrail for patient stability
- Large, 1" LCD Display
- Motion-sensing weighing technology
- Reweigh, LB/KG, Zero, On/Off Buttons. Low Battery Indicator and Auto-Off
- Heavy-duty steel base
- 19 1/8" x 19 1/8" platform with wheels for mobility
- 6 C-cell batteries & 120V AC adapter
- **3 Year Warranty**

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
1101KL	700 lb/310 kg	0.2 lb/0.1 kg	37 lb	19 1/8" x 19 1/8"	23 3/4" x 19 1/8" x 44 7/8"

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
1101KL	027045642422	1	UPC	11" x 43 1/4" x 19 3/4"	5.44	40 lb

Health o meter®
Professional