

## Digital Chair Scale

Our electronic chair scale is ideal for those who cannot stand on a conventional scale and who do not use a wheelchair. Durable and accurate, it also offers swing away arms and a fold up footrest to facilitate safe patient entry and exit.

### Features include:

- Motion-sensing weighing technology for improved accuracy and reliability
- Industrial-strength frame with sturdy wheels for mobility

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

[www.healthometermedical.com](http://www.healthometermedical.com)



Health o meter®  
Professional

# Electronic Chair Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)



595KL

- Motion-sensing weighing technology
- Large 1" LCD Display
- Reweigh, LB/KG, Zero, On/Off buttons. Low Battery Indicator and Auto-Off features
- 6 C-cell batteries (not included) & 120V AC adapter power supply (included)
- Industrial-strength steel frame
- Sanitary cushioned seat
- Movable arms for easy access
- Folding foot rest for patient comfort during weighing
- Heavy-duty wheels for mobility
- 3 Year Warranty

## Product Specifications

Model	Capacity	Increment	Product Wt.	Seat Dims	Product Dims
595KL	600 lb/270 kg	0.2 lb/0.1kg	79 lb	17 1/4"W x 17 3/4"D	23 1/4" x 35" x 40 1/8"H

## Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
595KL	027045642446	1	UPC	28" x 28" x 32"	4.15	85 lb

Health o meter®  
Professional