

# Pediatric and Neonatal

## ProPlus™ Neonatal Digital Pediatric Scale

By applying our leading-edge technology, Health o meter® has made weighing babies even more precise. Patient movement is compensated for and actual weights, not averages, are given. It's the ultimate in accuracy and reliability.



### Features include:

- Motion-sensing weighing technology compensates for infant movement to give actual—not average—weights
- Sensitive resolution (0.1 oz., 2g) for the neonatal environment

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

[www.healthometermedical.com](http://www.healthometermedical.com)

**Health o meter®**  
Professional

# ProPlus™ Neonatal Digital Pediatric Scale

Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)


Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)



2200KL

- Sensitive resolution (0.1 oz., 2 g) for the neonatal environment
- Motion-sensing weighing technology
- Reweigh, LB/KG, zero/tare and on/off buttons
- Auto-off and low-battery indicator features
- Easy-to-read 1" LCD display reads in pounds or kilograms with the touch of a button
- Removable tray
- Easy-to-clean sanitary acrylic tray with built-in measuring tape (up to 21 3/4" x 1/8" or 55 cm x 1/4 cm)
- Dual power source with rechargeable battery pack and 120V AC adapter (both included)
- 1 Year Limited Warranty

## Product Specifications

Model	Capacity	Increment	Product Wt.	Tray Dims	Product Dims
2200KL	0-12 lb/12 -30 lb 0-6 kg/6-15 kg	0.1 oz/0.2 oz 2 g/5 g	21 lb	12" x 22" x 4 1/4"	22" x 16" x 7 7/8"

## Packaging Specifications

Model	UPC	Case Pack	12 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
2200KL	027045642439	1	UPC	28 3/4" x 19" x 11"	3.48	25 lb

**Health o meter®**  
Professional