

Physician Digital Scales

For the ultimate in measurement accuracy, trust our physician digital scales. They have high weighing capacities, low profiles, large bases and user friendly interface.

Features include:

- 350° Swivel Display with easy-read 1" LCD



- 3 Styles: Eye-Level, Waist-High and Remote Display 2-Piece Platform
- Computer Interface

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



597KL
Eye-Level

Health o meter®
P r o f e s s i o n a l

Physician Digital Scales



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



- **BMI (easy scroll method)**
- **Hold/Release, Zero, KG/LB, Clear, Enter functions**
- **Auto-Zero, 5 minute auto-off (AC mode disabled), low battery indicator**
- **Wheels and Handle for maximum portability**
- **Silent Slide® height rod (24"-84" x 1/8" / 60cm – 213cm x 1mm) 597KL only**
- **Updatable Firmware**
- **RS232 Output Port**
- **6 C-cell batteries (not incl.) or 120V AC adapter (incl.)**
- **Metallic gray**

Additional Features for 752KL

- **Remote Display with 1" LCD and easy-touch keypad**
- **Fixture for wall-mount or table top display**
- **9' cord**

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
597KL	500 lb/227 kg	0.2 lb/0.1 kg	36 lb	14 1/4" x 14 1/4"	14 1/4" x 22" x 54 3/4"
599KL	500 lb/227 kg	0.2 lb/0.1 kg	28 lb	14 1/4" x 14 1/4"	14 1/4" x 22" x 44 1/8"
752KL	500 lb/227 kg	0.2 lb/0.1 kg	16 lb	14 1/4" x 14 1/4"	14 1/4" x 14 1/4" x 2 1/2"

Packaging Specifications

Model	UPC	Case Pack	Master Carton Dims	Case Cube	Package Wt.
597KL	022319408377	1	52 1/4" x 24 1/2" x 4 3/4"	3.52	38 lb
599KL	022319408360	1	40" x 21" x 4 3/4"	2.31	30 lb
752KL	022319408384	1	25 1/2" x 19 1/4" x 4 3/4"	1.35	18 lb

Health o meter®
Professional