



Pro Series Portable Digital Pediatric Scales

Health o meter® Professional has applied the accuracy of digital technology to pediatric scales that are also portable. With an easy-to-read display and factory calibration you're guaranteed precise weight readings even after the child is removed.

Features include:

- Load cell technology for improved accuracy and reliability
- Factory calibrated with in-the-field calibration
- Easy to read LCD display

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



549KL

Health o meter®
P r o f e s s i o n a l

Pro Series Portable Digital Pediatric Scale



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place
 Bridgeview, Illinois 60455
 USA Telephone: +1 (708) 598-9100
 Facsimile: +1 (708) 233-5472

Customer Service


US and Canada

Telephone: (800) 815-6615
 Facsimile: (800) 654-7330
HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100
 Facsimile: (708) 599-0150
HomProCS@pelstarllc.com



550KLS-01

- Load cell technology for improved accuracy and reliability
- Factory calibrated to assure accuracy; in-the-field calibration (requires 40 lb calibrated weight)
- Easy-to-read 0.5" LCD display reads in pounds or kilograms
- Two-minute hold feature for reading after child is removed.
- Easy-to-clean removable sanitary tray
- Dual power source with 9V battery (not included) or 120V AC adapter (available separately)
- Two minute auto-off prolongs battery life
- Convenient carrying case available
- **1 Year Limited Warranty**



549KL

- Load cell technology
- Factory calibrated to assure accuracy; in-the-field calibration (requires 10 kg & 20 kg calibrated weights)
- Easy-to-read 0.6" LCD display reads in pounds or kilograms with the touch of a button
- Easy-to-clean sanitary tray
- Powered by 1- 9V battery (included)
- Two minute auto-off prolongs battery life
- **1 Year Limited Warranty**

Product Specifications

| Model | Capacity | Increment | Product Wt. | Tray Dims | Product Dims |
|-----------|-------------|----------------|-------------|----------------------------|----------------------------|
| 550KLS-01 | 44 lb/20 kg | 1/2 oz/10 g | 8 lb | 21" x 9" x 3" | 21" x 11" x 4 1/2" |
| 549KL | 44 lb/20 kg | 1/2 oz/0.01 kg | 4 lb | 19 1/4" x 12 1/4" x 2 1/4" | 19 1/4" x 13 1/4" x 4 1/2" |

Packaging Specifications

| Model | UPC | Case Pk | I2 of 5 Code | Master Carton Dims | Case Cube | Pkg Wt. |
|-----------|--------------|---------|----------------|------------------------|-----------|---------|
| 550KLS-01 | 080006005247 | 1 | 20080006005241 | 22" x 13 1/4" x 7 1/2" | 1.27 | 10 lb |
| 549KL | 022319408308 | 1 | UPC | 21 1/2" x 15" x 5 1/4" | 0.98 | 5 lb |

