



## Pro Series Portable Digital Pediatric Scales

Health o meter® Professional has applied the accuracy of digital technology to pediatric scales that are also portable. With an easy-to-read display and factory calibration you're guaranteed precise weight readings even after the child is removed.

### Features include:

- Load cell technology for improved accuracy and reliability
- Factory calibrated with in-the-field calibration
- Easy to read LCD display

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

[www.healthometermedical.com](http://www.healthometermedical.com)



549KL

**Health o meter®**  
P r o f e s s i o n a l

# Pro Series Portable Digital Pediatric Scale



## Health o meter®

### Professional Products

#### Pelstar LLC

#### Corporate Headquarters

7400 West 100th Place  
 Bridgeview, Illinois 60455  
 USA Telephone: +1 (708) 598-9100  
 Facsimile: +1 (708) 233-5472

#### Customer Service

##### US and Canada

Telephone: (800) 815-6615  
 Facsimile: (800) 654-7330  
[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)

#### Customer Service

##### outside the US and Canada

Telephone: (708) 598-9100  
 Facsimile: (708) 599-0150  
[HomProCS@pelstarllc.com](mailto:HomProCS@pelstarllc.com)



#### 550KLS-01

- Load cell technology for improved accuracy and reliability
- Factory calibrated to assure accuracy; in-the-field calibration (requires 40 lb calibrated weight)
- Easy-to-read 0.5" LCD display reads in pounds or kilograms
- Two-minute hold feature for reading after child is removed.
- Easy-to-clean removable sanitary tray
- Dual power source with 9V battery (not included) or 120V AC adapter (available separately)
- Two minute auto-off prolongs battery life
- Convenient carrying case available
- **1 Year Limited Warranty**

#### 549KL

- Load cell technology
- Factory calibrated to assure accuracy; in-the-field calibration (requires 10 kg & 20 kg calibrated weights)
- Easy-to-read 0.6" LCD display reads in pounds or kilograms with the touch of a button
- Easy-to-clean sanitary tray
- Powered by 1- 9V battery (included)
- Two minute auto-off prolongs battery life
- **1 Year Limited Warranty**

### Product Specifications

Model	Capacity	Increment	Product Wt.	Tray Dims	Product Dims
550KLS-01	44 lb/20 kg	1/2 oz/10 g	8 lb	21" x 9" x 3"	21" x 11" x 4 1/2"
549KL	44 lb/20 kg	1/2 oz/0.01 kg	4 lb	19 1/4" x 12 1/4" x 2 1/4"	19 1/4" x 13 1/4" x 4 1/2"

### Packaging Specifications

Model	UPC	Case Pk	I2 of 5 Code	Master Carton Dims	Case Cube	Pkg Wt.
550KLS-01	080006005247	1	20080006005241	22" x 13 1/4" x 7 1/2"	1.27	10 lb
549KL	022319408308	1	UPC	21 1/2" x 15" x 5 1/4"	0.98	5 lb

