

ProPlus® Digital Wheelchair Ramp Scales

Compensating for movement is critical to weighing patients in wheelchairs. That's why this scale has the exclusive technology to do just that. Extremely durable, with state-of-the-art technology, these wheelchair scales are a "must-have" for any long-term and acute-care facility.

Features include:

- 180° swivel head for side read
- Customizable readout
- Keypad entry.
- Large platforms can be used with chairs, wheelchairs or a stand on scale
- Computer interface with 270 patient memory

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



2500KL

Health o meter®
Professional

ProPlus® Wheelchair Ramp Scales



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



2500KL



2600KL

- 180° Swivel Head for Side Read
- Motion-sensing weighing technology
- Customizable 1" backlight LCD display
- Displays fractions, lb/oz, decimals in standard; decimals in metric
- User-friendly menu screen
- Platform with wheels for mobility, 2600KL folds up for easy storage
- Keypad entry, BMI, Manual Tare, Zero, Auto-Off, Auto-Hold, Reweigh, Hold/Release
- USB and Serial Ports (2) for computer interface with 270 patient memory
- Updatable firmware
- 6 D-cell batteries (6 month battery life) & 120V AC adapter
- Ramp attaches to either side for left/right access. Second ramp available (2600KL only)
- 3 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
2500KL	800 lb 360 kg	4 oz/0.2 lb/1/4 lb 0.1 kg	71 lb	28 1/2" x 29 3/8" (10" ramp)	39 1/2" x 37 7/8" x 46 1/2" (ramp down)
2600KL	800 lb 360 kg	4 oz/0.2 lb/1/4 lb 0.1 kg	105 lb	32" x 37" (10" ramp)	42" x 45" x 46 1/2" (ramp down)
2600RAMP				10" x 32"	

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
2500KL	027045642477	1	UPC	45 1/2" x 34" x 9 7/8"	8.82	82 lb
2600KL	027045642484	1	UPC	13 3/4" x 35 1/2" x 40"	11.22	10.5
2600RAMP		1	UPC			

Health o meter®
Professional