

Balance Beam Wheelchair Scales

Health o meter's strength and accuracy is fully evident in this unique specialized scale. With its large steel platform for oversized wheelchairs, poise bar lock-in and ramp attachments on either side, taking measurements is easy.

Features include:

- Rotating poise bars for weighing in pounds and kilograms. Poise bar lock-in for error free operation
- Permanent, engraved readings on poise bar for long lasting readability

For a complete set of features, specifications and packaging information, see reverse side.

Health o meter® Professional is the industry leader in medical scales providing advanced innovations and the highest level of quality and accuracy since 1919.

www.healthometermedical.com



410KLS-01

Health o meter®
P r o f e s s i o n a l



Health o meter®

Professional Products

Pelstar LLC

Corporate Headquarters

7400 West 100th Place

Bridgeview, Illinois 60455

USA Telephone: +1 (708) 598-9100

Facsimile: +1 (708) 233-5472

Customer Service

US and Canada

Telephone: (800) 815-6615

Facsimile: (800) 654-7330

HomProCS@pelstarllc.com

Customer Service

outside the US and Canada

Telephone: (708) 598-9100

Facsimile: (708) 599-0150

HomProCS@pelstarllc.com



421KLS-01
(ramp only
available as
420S-01)



410KLS-01

- Consists of a balance beam scale with large 30" x 30" steel platform (410KLS-01) to accommodate oversized wheelchairs
- Ramp attaches to either side for left or right access (410KLS-01)
- Rotating poise bars for weighing in pounds and kilograms. Poise bar lock-in for error-free operation
- Permanent, engraved readings on poise bar for long-lasting readability
- Heavy-duty steel base for improved strength and durability
- Factory calibrated to ensure accuracy
- Includes two counterweights to increase capacity to 500 lb or 200 kg
- 2 Year Limited Warranty

Product Specifications

Model	Capacity	Increment	Product Wt.	Platform Dims	Product Dims
410KLS-01	500 lb/200 kg	1/4 lb/100 g	175 lb	30" x 30"	30" x 30" x 57 1/4"
421KLS-01	500 lb/200 kg	1/4 lb/100 g	75 lb	26" x 28 3/8"	26 1/4" x 33 3/4" x 57 1/4"
420S-01	Ramp only				

Packaging Specifications

Model	UPC	Case Pack	I2 of 5 Code	Master Carton Dims	Case Cube	Package Wt.
410KLS-01	080006008668	1	UPC	60" x 34 1/4" x 32"	38.06	185 lbs
421KLS-01	080006011200	1	UPC	59 1/4" x 6 3/4" x 20 1/2" (Scale)	4.74	38 lbs
				34" x 28" x 5" (Ramp)	2.75	47 lbs
420S-01	080006008675	1	20020006008679	34" x 28" x 5"	2.75	47 lbs

Health o meter®
Professional