

High Sensitivity Probes

Description:

Based on over 20 years experience in this field, the latest generation of the world renown Dopplex® handheld Dopplers offer even greater performance, quality and value.

Probe improvements include:

- New design with 50% greater sensitivity
- New EZ8 wide beam probe for easy vessel detection
- All fixed obstetric probes are now waterproof

EZ 8 Probe: The new EZ8 High Sensitivity Doppler probe incorporates wide beam technology to allow easy location of the vessel. It is also easier to maintain vessel contact during inflation and deflation.




VP4HS: A 4MHz High Sensitivity Doppler probe for detection of deep lying vessels.




VP5HS: A 5MHz High Sensitivity Doppler probe for edematous limbs and deep lying vessels. The ideal probe as an adjunct to the EZ8 for ABI measurements.



VP8HS: An 8MHz High Sensitivity Doppler probe for easy detection of peripheral vessels and calcified arteries.



VP10HS: A 10MHz High Sensitivity Doppler probe for detecting smaller vessels in superficial applications.



OP2HS: A 2MHz High Sensitivity Doppler probe offers enhanced performance in 2nd & 3rd trimester pregnancy



OP3HS: A 3MHz High Sensitivity Doppler probe with increased sensitivity to low velocity signals. For use in early pregnancy, as early as 8-10 weeks gestation.



New High Sensitivity Vascular Probes



List Price 2003

\$316.00 Standard Probes

\$346.00 EZ8 Wide Beam Probe