

LESS IS MORE!

- **LESS TO DO**

Fewer steps to accomplish the same procedures

- **LESS TO BUY**

Fewer accessories needed for additional functionality

- **LESS TO PAY**

Lower initial costs than competitors, with lower priced accessories

MORE USER FRIENDLY



Hand Control Jack For Low And High Power

Return Electrode Jack

Bipolar Jack

Foot Control Jack

tTJT -0.0001 Tc0.0001 Tw (anps (A80)16.50w).9(1) t.s liquid

Change Power Advantage – Less To Do...

The Aaron 900™ allows the doctor to change the power setting with one action, using the rotary power-control dial or the up/down buttons on the autoclavable handpiece.

Here are the steps necessary just to go from Low to High Power on the Hyfreator® 2000.

(reverse steps to go from high to low power)



Step 1

Unplug the handpiece from the Lo power receptacle



Step 2

Move the power selector switch from Lo to Hi




Step 3

Plug the handpiece into the Hi power receptacle



Step 4

Adjust the rotary dial or up/down buttons on the handpiece



Technology Advantages – Less To Worry About...



The Bovie 3-Button Handpiece – Better By Design, Inside And Out

The basic function of a 3-button handpiece is to allow the user to control power up/down and to activate the high frequency desiccator. The Aaron 900 has a technology advantage right in the handpiece: a superbly engineered collet that's even gold plated for superior conductivity. The collet is attached to our circuit board and the entire handpiece is encapsulated and formed in a high tech polymer that combines comfort for you with a sure grip. The buttons on the handpiece have been designed for superior tactile feel and comfort. Finally, it's autoclavable so you can have sterility when you need it, and an easy to clean and durable handpiece all the way around.

Digital Error Detection

Digital error detection means unsurpassed safety for you and your patient. The Aaron 900 continuously monitors every aspect of the unit's output. At the sign of any problem the machine instantly disables the output and displays the appropriate error code on the display.



SPECS

- Output Power 0-10 Watts adjustable in 1/10 Watt increments
10-30 Watts adjustable in 1 Watt increments
- Line Voltage A900 (A800EU-110) – 120 VAC ± 10%
A900-220 (A800EU-220) – 220 VAC ± 10%
- Output Frequency 550 kHz
- Line Frequency 50-60 Hz
- Dimensions Width = 22.2 cm (9") • Height = 15.9 cm (6.25")
Depth = 11.4 cm (4.5") • Weight = < 2.3 kg (< 5 lbs.)
- Warranty Two Years (unit only)

Designed For Years of Trouble-Free Use

The Aaron 900 is designed and engineered to be the most reliable and durable desiccator available today. The unit's life is prolonged by incorporating the latest in chip technology and reducing the number of internal parts, thereby generating less internal heat. Post-procedure and long-term cleanup are simplified by its non-absorbing, nonstaining, high-impact plastic case.

Package



Unit, power cord, handpiece (A901), 50 non-sterile disposable sharp dermal tips (A804), 50 non-sterile disposable blunt dermal tips (A806), wall mount kit (A837), and user's guide.

A Lower Cost Alternative

If cost is a major factor, consider the Aaron 800-EU. You get the same great technology, just without the 3-button handpiece and with fewer included electrodes. It costs around \$100 less.

Package includes unit, power cord, handpiece, one non-sterile reusable electrode, wall mount kit, and user's guide.



DISTRIBUTED BY:

Bovie® **AARON**
MEDICAL®

7100 30th Ave. N. • St. Petersburg, FL 33710-2902
U.S. Phone 1-800-537-2790 • Fax 1-800-323-1640
Int'l. Phone +1-727-384-2323 • Fax +1-727-344-3876
www.aaronmed.com • sales@aaronmed.com