

EPIC SERIES

Clinton's Epic Series combines comfort, functional features and a large storage area to bring you one of the most versatile tables in our line. Our Epic Series table comes in your choice of 14 colors of upholstery and color coordinating doors on sturdy, easy-clean cabinet base.

FEATURES

- Plush vinyl upholstery over 2" firm padding
- Full 27" width for patient comfort
- Air-spring backrest
- Concealed paper dispenser
- Nylon coated, two position stirrups
- Pull-out leg extension
- Pull-out basin shelf
- Easy care all laminate base
- Color coordinated doors and top
- Chip resistant PVC edges
- Two large storage cabinets, each with an adjustable shelf
- Enamel coated, welded steel top frame
- Tables can be made right or left-handed
- 350 lbs. capacity
- 3-year limited warranty



Adjustable Stirrups

Air Lift Back



Photo above shows cabinet's easy clean laminate interior and the hide-away stirrups in the extended position.

**ORDERING INFORMATION
FOR MODEL # 8840**

1. Select base color (page 37 for legs and frame)
2. Choose a door color from the Color Selector on page 37
3. Pick a vinyl color from the Color Selector on page 37

	Length	Height	Width
8840	56"/70" extended	30.5"	27"

- Full cabinet model table with two doors and adjustable shelves
- Can be adjusted for either right or left-handed installation
- Two position stirrups, air-spring backrest, padded leg extension and basin shelf are standard
- If damage occurs, all components are easily replaceable in the field
- Modular construction allows for easy delivery by one person