


ETA STYLE LINE SERIES

Contemporary design and easy clean surfaces highlight the ETA (Easy-To-Assemble) Style Line Series of tables. Its flat panel components use bolt-through, pre-threaded assembly hardware for an easy and secure assembly. Plus, Clinton's Style Line Series is available in 10 coordinated laminate and vinyl colors to compliment any decor.

FEATURES

- Save money on shipping. Clinton ETA tables ship flat by UPS® or FedEx® Ground (1 to 3 units)
- Single person delivery and assembly
- No inside delivery charges (UPS® or FedEx® Ground)
- 350 lbs. load capacity
- Assemble with only basic tools
- Large orders that are less expensive to ship via common carrier will have the tops pre-attached (4 or more units)

Plywood frame with steel reinforced stress points.



	Length	Height	Width
9001-27	72"	31"	27"
9001-30	72"	31"	30"

Straight Line Treatment Table

- Features all laminate surfaces and adjustable backrest
- Quick Ship in Gray and Natural only.



	Length	Height	Width
9004-27	72"	31"	27"
9004-30	72"	31"	30"

Straight Line Treatment Table

- Features shelf with storage compartments and adjustable backrest.
- Quick Ship in Gray and Natural only.



	Length	Height	Width
9013-27	72"	31"	27"
9013-30	72"	31"	30"

Straight Line Treatment Table

- Features two drawers, shelf and adjustable backrest
- Quick Ship in Gray and Natural only.



	Length	Height	Width
9020-27	72"	31"	27"
9020-30	72"	31"	30"

Straight Line Treatment Table

- Features full length storage shelf and adjustable backrest
- Quick Ship in Gray and Natural only.
- Color coordinated shelf must be specified

