

	Length	Height	Width
<b>1702-27</b>	48"	36"	27"
<b>1702-30</b>	48"	36"	30"

### Taping Table

- Heavy duty taping table with solid wood legs, H-braces and laminated full shelf
- 2" extra firm foam padding
- 350 lbs. capacity

	Length	Height	Width
<b>1701-27</b>	48"	36"	27"
<b>1701-30</b>	48"	36"	30"

### Taping Table (not shown)

- Same as 1702 without shelf

	Length	Height	Width
<b>1764</b>	78"	31"	30"

### Double Leg Lift Table

- Heavy duty table with panel legs, 2 level laminated storage compartments and shelf
- Each leg rest can be raised and lowered independently in 1 1/4" increments
- 2" extra firm foam padding
- 400 lbs. capacity

	Length	Height	Width
<b>1703-27</b>	48"	36"	27"
<b>1703-30</b>	48"	36"	30"

### Taping Table

- Heavy duty taping table with solid wood legs, H-braces and laminated 2 level storage compartments and shelf
- 2" extra firm foam padding
- 350 lbs. capacity

	Length	Height	Width
<b>1763</b>	78"	31"	30"

### Knee Gatch Table

- Heavy duty table with panel legs, 2 level laminated storage compartments and open shelf
- 2" foam padding
- Adjustable backrest
- Adjustable foot section
- 350 lbs. capacity

	Length	Height	Depth
<b>1772</b>	76"	36"	30"

**2 Station Taping Table**

- All natural oak plywood construction with a durable finish
- 2" extra firm foam padding
- 2 padded 27" x 30" tops
- 3 storage sections, each with an adjustable shelf
- Self-closing doors with fully adjustable Euro-style hinges
- 1 large laminate center drawer
- 2 open storage shelves
- Edgebanding on all edges
- **\* Custom models available**
- 350 lbs. load capacity each padded seat

	Length	Height	Width
<b>1780—table</b>	78"	36"	30"
	Depth	Height	Width
back section	30"	36"	30"

**3 Station Taping Table/Treatment Table Combo**

- All natural oak plywood construction with a durable finish
- 2" extra firm foam padding
- Dual purpose table—use flat as a treatment table or as 3 individual taping stations
- 3 open storage sections, each with 2 adjustable shelves
- 1 enclosed cabinet with 2 doors and 1 adjustable shelf
- Edgebanding on all edges
- 350 lbs. load capacity

Back View

	Length	Seat Height	Depth
<b>1790</b>	76"	36"	51"

**Multi-Position Training Table**

- Same features as model 1772 with 30" high stationary backrest and 51" depth