

STYLE LINE FULL CABINET MODELS


Contemporary design and easy clean surfaces highlight the Style Line Series of tables. Great looking and durable, Clinton's Style Line Series is available in 10 coordinated laminate and vinyl colors to complement any decor.

FEATURES

- All cabinet models ship fully assembled
- Large selection of models to fit your budget
- Specialty models to fit your needs
- 350 lbs. load capacity
- Color coordinated door laminate and vinyl available
- All models ship by common carrier only

Hardwood plywood frame with steel reinforced stress points.

Heavy-duty, all steel adjustable backrest support system



	Length	Height	Width
9070-27	72"	31"	27"
9070-30	72"	31"	30"

Cabinet Style Laminate Treatment Table

- Features 2 sliding doors front and back plus adjustable backrest
- Quick Ship in Gray and Natural only. **Quick Ship**

	Length	Height	Width
9044-27	72"	31"	27"
9044-30	72"	31"	30"

Cabinet Style Laminate Treatment Table

- Features an adjustable backrest, 2 storage compartments with sliding doors and an adjustable shelf in each compartment

	Length	Height	Width
9105-30	72"	31"	30"

Cabinet Style Laminate Treatment Table

- Features 3 drawers, adjustable backrest and 1 storage compartment with 1 adjustable shelf

	Length	Height	Width
9074-24	72"	31"	24"
9074-27	72"	31"	27"
9074-30	72"	31"	30"

Cabinet Style Laminate Treatment Table

- Features four doors, adjustable backrest and 2 storage compartments each with an adjustable shelf