

STYLE LINE SPECIALITY TABLES

	Length	Height	Width	Foot Rest Width
9110-24	72" extended/56" folded	31"	24"	21"
9110-27	72" extended/56" folded	31"	27"	23"

Space Saver Table

- Features 1 adjustable shelf, 2 doors, adjustable footrest and adjustable backrest
- 350 lbs. load capacity

	Length	Height	Width
9107-38	78"	31"	30"

Cabinet Style Laminate Treatment Table

- Features 3 drawers, 2 storage compartments with 2 adjustable shelves and adjustable backrest
- 350 lbs. load capacity

	Length	Height	Width
9133-27	72"	31"	27"
9133-30	72"	31"	30"

Laminate Treatment Table with Stool

- Features 3 drawers, storage shelf, adjustable backrest and matching stool that stores neatly under the shelf
- Stool is 18" wide
- 350 lbs. load capacity

	Length	Height	Width
9102-30	72"	31"	30"

Straight Line Treatment Table

- Features 14" wide double pull-out shelves (accessed from one side only) and a full length shelf
- Shipped KD By common carrier only
- 350 lbs. load capacity

	Length	Height	Width
9104-38	78"	31"	30"

Special Procedure Table with 2 Backrests

- Features cabinet style laminate table with 2 enclosed storage compartments with adjustable shelves, 22" wide pull-out shelf with locking length positions and 2 adjustable backrests
- 350 lbs. load capacity

	Length	Height	Width
9106-38	78"	31"	30"

Laminate Cast Table

- Features adjustable backrest and 3 large storage sections with adjustable shelves
- 350 lbs. load capacity