

## REPLACEMENT PARTS

### 158 Replacement Batteries

*Our complete line of replacement batteries will keep your Welch Allyn medical products working for years to come.*

### 159 Replacement Lamps

*We offer replacement lamps to keep your Welch Allyn instruments working and dependable.*




























### 167 Replacement Components for Products

*A comprehensive listing of disposable replacement items — all in one easy-to-find place.*



















# REPLACEMENT LAMPS

Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments
00200		<p><b>Otosopes</b></p>  20100  21600	<p><b>Otosopes PenLite</b></p>  235  77900
		<p><b>Anoscopes</b></p>  38800  380 Series  385 Series  395 Series	
		<p><b>Illuminators, Transilluminators</b></p>  26030  41000	<p><b>Illuminators, Transilluminators</b></p>  26000  42700  40510
00600			<p><b>Sigmoidoscopes</b></p> <p>30800</p>
00900		<p><b>Retinoscopes</b></p>	<p><b>Ophthalmoscopes Retinoscopes</b></p>  10600  121  110  17600  17610
01300		<p><b>Tongue Blade Holders</b></p>  28000	<p><b>PenLites</b></p>  76000  77700
01400			<p><b>PenLite</b></p>  77800






























Diagnostic Instruments  
 Diagnostic — Vital Signs  
 ECG/Stress/Holter/Spirometry  
 Eye & Ear Care  
 Laryngoscopes  
 Lighting  
 Monitoring Devices  
 Women's Health/Endoscopy  
 Replacement Parts  
 Index/Reference

# REPLACEMENT LAMPS

Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments
02000			<b>Headlights</b>  45003
02200			<b>Otosopes</b>  22900
02500		<b>Headlights</b>  46003	
02600		<b>Headlights</b>  49003	
03000		<b>Ophthalmoscopes</b>  11710 <b>StrabismoScope™</b>  12400 <b>EpiScope™</b>  47300	<b>Ophthalmoscopes</b>  11600 11605 11610 <b>Retinoscopes</b>  18000










Diagnostic Instruments  
 Diagnostic — Vital Signs  
 ECG/Stress/Holter/Spirometry  
 Eye & Ear Care  
 Laryngoscopes  
 Lighting  
 Monitoring Devices  
 Women's Health/Endoscopy  
 Replacement Parts  
 Index/Reference

# REPLACEMENT LAMPS

Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments	
03100		<p><b>Otosopes</b></p>  20000  25020  21700  3.5v Handle Adapter 73500	<p><b>Otosopes</b></p>  25000  25200	
		<p> 20200               20202               Tongue Blade Holder 28100</p> <p><b>Illuminators, Transilluminators</b></p>  27000  27050  26530  41100  43300	26500	
03300 (02100)			<p><b>Ophthalmoscopes</b></p>  11511  11500	
03400 (03800) (02700)		<p><b>Otosopes</b></p>  21110  24031	<p><b>Otosopes</b></p>  24011  24020	
		<p><b>2.5V Handle Adapter</b></p>  73550	<p><b>Illuminators</b></p>  27200  27250	<p><b>Otosopes</b></p>  24000
03600		<p><b>Disposable Laryngoscope Illuminator</b></p>  60500		













Diagnostic Instruments  
 Diagnostic — Vital Signs  
 ECG/Stress/Holter/Spirometry  
 Eye & Ear Care  
 Laryngoscopes  
 Lighting  
 Monitoring Devices  
 Women's Health/Endoscopy  
 Replacement Parts  
 Index/Reference

# REPLACEMENT LAMPS



Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments
03700			<b>Retinoscopes</b>  18100
03800		<b>PanOptic™ Ophthalmoscope</b>  11810      11820      11800	
03900		<b>Ophthalmoscopes</b>  12810	<b>Ophthalmoscopes</b>  11110
04050		<b>LX-150 Light Source</b>  45250	
04100			<b>Examination Lights</b>

Diagnostic Instruments  
 Diagnostic — Vital Signs  
 ECG/Stress/Holter/Spirometry  
 Eye & Ear Care  
 Laryngoscopes  
 Lighting  
 Monitoring Devices  
 Women's Health/Endoscopy  
 Replacement Parts  
 Index/Reference















# REPLACEMENT LAMPS

Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments
04200	 (shown at 75% of actual size)	<b>Examination Lights</b> 	48600 48610 48625 48635
04400 (04300)		<b>Ophthalmoscopes</b> 	<b>Ophthalmoscopes</b> 
04500			<b>Retinoscopes</b> 
04600			<b>Retinoscopes</b> 
04700		<b>Standard Laryngoscopes</b> 	 684 Series Sizes 0-1                      694 Series Sizes 1-2

Diagnostic Instruments  
 Diagnostic — Vital Signs  
 ECG/Stress/Holter/Spirometry  
 Eye & Ear Care  
 Laryngoscopes  
 Lighting  
 Monitoring Devices  
 Women's Health/Endoscopy  
 Replacement Parts  
 Index/Reference









# REPLACEMENT LAMPS

Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments	
06300			<b>Examination Lights</b> 	
06400	 (shown at 75% of actual size)	<b>Examination Lights</b>  44200		
07200		<b>Retinoscope</b>  17200		
07600		<b>PenLite</b>  76600		
07800 (00800) (04000) (01900)		<b>Vaginal Specula Illuminator</b>  78000	<b>Sigmoidoscope, Anoscope Light Handle</b>  73210	<b>Sigmoidoscope, Anoscope Light Handle</b>  73200
08100			<b>Examination Lights</b>  48100	

Diagnostic Instruments  
 Diagnostic — Vital Signs  
 ECG/Stress/Holter/Spirometry  
 Eye & Ear Care  
 Laryngoscopes  
 Lighting  
 Monitoring Devices  
 Women's Health/Endoscopy  
 Replacement Parts  
 Index/Reference

# REPLACEMENT LAMPS

Lamp Number	Lamp-Actual Size	Current Instruments	Obsolete Instruments
08200		<b>Retinoscopes</b>  18200	
08300		<b>Retinoscopes</b>  18300	
09500	 (shown at 25% of actual size)		
09600			

All Lamps (except 02000, 02500, 02600, 04050, 04060, 04100, 04200, 06100, 06400, 08100, 09500, 09600 and 09800) Sold to Distributors in Box Quantity only — 6 Lamps/Box.

### Welch Allyn Lamps Policy

Welch Allyn instruments are completely engineered from the source of illumination to the instrument's usage. A change in this sequence reduces product effectiveness. Lamp manufacture is a very critical part of this system, therefore, we recommend that you use only genuine Welch Allyn lamps.

### Welch Allyn Batteries Policy

A Welch Allyn nickel cadmium battery will be replaced on a pro rata basis should it fail before the date printed on the battery (when used in Welch Allyn handles only). This date represents two years after manufacture of the battery.

## Replacement Components Not Included in Charts

*Common replacement components in one place.*

If you don't find what you need here, check out the replacement parts catalog online at [www.welchallyn.com/technical](http://www.welchallyn.com/technical) service.

### Accessory Numbers

#### Disposable Probe Covers for Thermometers

- 05031-110: Disposable Probe Covers (10,000)
- 05051-001: Disposable Probe Tips (1,000)
- 05051-005: Disposable Probe Tips (5,000)

#### Disposable Otoscope Specula

- 52133: 521 Series KleenSpec® Disposable Otoscope Specula, 3 mm
- 52134: 521 Series KleenSpec® Disposable Otoscope Specula, 4 mm
- 52135: 521 Series KleenSpec® Disposable Otoscope Specula, 5 mm
- 52432: 524 Series KleenSpec® Disposable Otoscope Specula, 2.5 mm
- 52434: 524 Series KleenSpec® Disposable Otoscope Specula, 4 mm
- 52450: KleenScoop™ Cerumen Removal Specula, 10 Boxes of 100 each
- 52450-B: KleenScoop™ Cerumen Removal Specula, Box of 100

#### Disposable Sheaths

- 52640: Disposable Sheaths for Exam Light (25/box)

#### Forms

- 52300: Audiogram form (1 pad of 50, 10 pads per case)
- 55230: AudioScope® Recording Forms (100/pkg, 10 pkg/box)

#### Disposable Flow Tubes

- 26001-0000: Disposable Flow Tubes, SP-150/SP-160 (10 each)
- 26001-0001: Disposable Flow Tubes, SP-150/SP-160 (100 each)

#### Disposable Nose Clips

- 56130: Disposable Nose Clips

#### Eartips

- 26008: Eartips for TM 262™, 8 mm, 1 box of 25
- 26011: Eartips for TM 262™, 11 mm, 1 box of 25
- 26013: Eartips for TM 262™, 13 mm, 1 box of 25
- 26015: Eartips for TM 262™, 15 mm, 1 box of 25
- 26017: Eartips for TM 262™, 17 mm, 1 box of 25
- 26019: Eartips for TM 262™, 19 mm, 1 box of 25

#### Mouthpieces

- 56120: Pediatric Disposable PneumoCheck™ Mouthpiece
- 56121: Adult Disposable PneumoCheck™ Mouthpiece
- 26002-0000: Cardboard Mouthpieces, SP-100/SP-110 (100/box)
- 26003-0000: Cardboard Mouthpieces, SP-20 (100/box)

#### Paper

- 52600: TM 262™ Printer Paper
- 53600: Printer Paper (4 rolls/box, 6 box/case)
- 56100: Thermal Paper (5 rolls)
- 7052-25: Printer Paper, 6 rolls
- 6200-40: Printer Paper for Atlas Monitor (box of 25)
- 94001-0000: Paper: SP-1, 3.5" x 3.5" (280 sheets/pack; 10 packs/case)
- 94002-0000: Paper: 8.5" x 11" (215 sheets/pack; 8 packs/case)

- 94003-0000: Paper: AT-3, 3.5" x 2.75" (175 sheets/pack; 12 packs/case)
- 94004-0000: Paper: AT-4, 3.5" x 3.5" (255 sheets/pack; 10 packs/case)
- 94005-0000: Paper: AT-5, 5.3" x 5.5" (140 sheets/pack; 12 packs/case)
- 94006-0000: Paper: AT-6, 5.3" x 5.5" (350 sheets/pack; 6 packs/case)
- 94008-0000: Paper: CS-6/12, 8.5" x 5.5" (350 sheets/pack; 6 packs/case)
- 94009-0000: Paper: AT-60, CS-100, 8.5" x 5.5" (250 sheets/pack; 4 packs/case)
- 94010-0000: Paper: AT-10, 8.5" x 5.5" (150 sheets/pack; 10 packs/case)
- 94011-0000: Paper: CM-8 (1 pack)

#### **Cuff Accessories**

- 5086-01: Premium Medium Inflation Bulb (Black)
- 5086-02: Rubber Bulb
- 5088-01: Air Release Valve and Bulb
- 5086-03H: Large Inflation Bulb, 6 Pack
- 5087-01: Air Release Valve
- 5091-13: Inflation System Basket (Small, Stainless Steel)
- 5090-65: Sphygmomanometer Repair Kit
- 5091-47: Inflation System Basket (Large, Black)
- 7670-07: Inflation System Basket, Large, White

#### **Batteries (not included in chart)**

- 72250: Rechargeable Battery for LumiView Portable Power Source
- 72500: 7.2v Rechargeable Battery (2)
- 72700: 7.2v Rechargeable Battery for Original MicroTymp®
- 72900: 2.4v Nickel-Cadmium Rechargeable Battery for MicroTymp® 2
- 6200-41: Replacement Battery for Atlas Monitor
- 72320: Rechargeable Battery Pack Only
- 72329: Complete Battery Pack Assembly with Rechargeable Batteries