

Recommended Inpatient Care Set-Ups

ROUTINE PHYSICAL ASSESSMENT PRODUCTS

Product Description	Part Number	Page
Universal Desk Charger Set	#71811-PSM	5
3.5v Diagnostic Set	#97200-C	6
Welch Allyn Disposable & Durable One-Piece Cuffs	(various)	9, 31
Professional Stethoscope	#5079-135	9
Uniscope—Disposable Stethoscope	#17462	9
767 Wall Aneroid	#7670-01	10
767 Mobile Aneroid	#7670-04	10
Spot Vital Signs®	#42NTB-E1	12
Welch Allyn Vital Signs Monitor 300 Series	#53NTP-E1	13
SureTemp® Plus 692 Electronic Thermometer	#01692-200	14
Braun ThermoScan® Pro3000 Thermometer	#03000-200	15
LS-150 Halogen Exam Light	#44500	16

PATIENT MONITORING AND DEFIBRILLATION PRODUCTS

Product Description	Part Number	Page
Propaq® CS Multiparameter Patient Monitor	#01031659	19
Acuity® and Acuity® LT Central Monitoring Stations	(custom)	20, 21
Micropaq® Patient-Worn Monitor	#01080006	21
PIC 30™ Monitor/Defibrillator	#973039C00	23
Atlas™ Monitor	#622NP-E1	24

SPECIALTY PROCEDURE PRODUCTS

Product Description	Part Number	Page
KleenSpec® Disposable Vaginal Specula and Light Source	#58001S, 78810	26
Fiber Optic Laryngoscope Handles & Blades	(various)	27
Exam Light III™ Fiber Optic Light	#48700	28
LS-200™ Halogen Procedure Light	#44200	28
Solid State Portable Headlight	#46070	28
Solid State Procedure Headlight	#49020	28



INPATIENT CARE

4341 State Street Road, P.O. Box 220, Skaneateles Falls, NY 13153-0220 USA
(P) 800-535-6663 (F) 315-685-3361 www.welchallyn.com

WelchAllyn®

Advancing Frontline Care™

©2003, Welch Allyn, Inc. Printed in USA SM2438 Rev B

Nellcor® is a registered trademark of Nellcor Puritan Bennett Inc. ThermoScan® is a registered trademark of Braun GmbH.

WelchAllyn®

Advancing Frontline Care™

Welch Allyn. The Point of Difference at the Point of Care.

It's the point where you and the patient meet. The point of contact, where you have to evaluate, monitor, and provide treatment—quickly and accurately, without losing sight of the patient as a person.

The very point, in fact, where we at Welch Allyn focus our skill. As the first name in medical diagnostic instruments, we understand the challenges you face in caring for your patients. The value of time, the importance of accuracy, simplicity and mobility, the ability to deliver the best care under demanding conditions—we've been there, for 88 years.

The result is innovation which puts state-of-the-art technology at your service, with practical, efficient, reliable, affordable instruments that give you the very best tools for a full spectrum of hospital inpatient applications.

It's the point where what we do best enhances what you do best.

The point of care.

INPATIENT CARE

- 1 Acuity® LT Central Monitoring Station Pg. 20
- 2 Universal Desk Charger Set #71641 Pg. 5
- 3 PanOptic™ Ophthalmoscope Head #11820 Pg. 7
- 4 SureTemp® Plus 692 Electronic Thermometer #01692-200 Pg. 14
- 5 Braun ThermoScan® Pro3000 #03000-200 Pg. 15
- 6 Fiber Optic Laryngoscope Handles & Blades Pg. 27
- 7 PIC 30™ Monitor/Defibrillator #973039CPO Pg. 23
- 8 AED 10™ #970300 Pg. 22
- 9 Tycos® 767-Series Mobile Aneroid #7670-04 Pg.10
- 10 LS-150 Halogen Exam Light #44500 Pg. 28
- 11 Spot Vital Signs® #42NTB-E1 Pg. 12
- 12 Professional Stethoscope #5079-135 Pg. 9
- 13 Vital Signs Monitor 300 Series #53NTP-E1 Pg. 13
- 14 Durable One-Piece Blood Pressure Cuffs Pg. 9
- 15 Propaq® CS Monitor, 242 Series #1031633 Pg. 19
- 16 Micropaq® #01080001 Pg. 21
- 17 Disposable Blood Pressure Cuffs Pg. 9

2 For more information on the products in this catalog, visit our website at www.welchallyn.com



ROUTINE PHYSICAL ASSESSMENT

“

”



3.5v Diagnostic Set

Keep your diagnostic instruments together and easily portable, with our most popular handle and instruments in one easy-to-order set.

- Choice of the PanOptic™ Ophthalmoscope or 3.5v Coaxial Ophthalmoscope (see facing page for details)
- Diagnostic Otoscope with Halogen HPX™ lamp provides 30% more light output
- Convertible Ni-Cad handle uses rechargeable battery (included) or C-cell batteries
- Comes complete with a zippered, hard case for secure storage and easy portability

#97200-C 3.5v Diagnostic Set with Coaxial Ophthalmoscope (#11720), Diagnostic Otoscope (#20000) and convertible Ni-Cad handle (#71000-C), includes specula assortment and hard zipper case

#97200-PC 3.5v Diagnostic Set with PanOptic™ Ophthalmoscope (#11820), Diagnostic Otoscope (#20000) and convertible Ni-Cad handle (#71000-C), includes specula assortment and hard zipper case

#97200-C


zipper case

767 Integrated Diagnostic Systems

- Choice of the revolutionary PanOptic™ Ophthalmoscope (see facing page) or a 3.5v Coaxial Ophthalmoscope
- Diagnostic Otoscope with Halogen HPX™ lamp provides 30% more light output
- Handle-based rheostats for quick light intensity adjustment
- 767 Wall Aneroid with easy-to-read dial face, certified accuracy, and lifetime calibration program
- KleenSpec® Plus Disposable Specula Dispenser, pre-filled with ear specula
- Optional locking collar available to prevent theft (#76700)
- Optional—767 SureTemp® Electronic Thermometer for highly accurate oral readings in just four seconds



All 767 Integrated Diagnostic Systems include Wall Transformer (#76710), Diagnostic Otoscope (#25020), 767 Wall Aneroid (#7670-01) and Disposable Specula Dispenser (#52401).

#76791-2P PanOptic™ Ophthalmoscope (#11820), 767 SureTemp® Thermometer (#76751)	#76792-2P PanOptic™ Ophthalmoscope (#11820), Thermometer not included	#76793 Coaxial Ophthalmoscope (#11720), 767 SureTemp® Thermometer (#76751)	#76794 Coaxial Ophthalmoscope (#11720), Thermometer not included
---	---	--	--



“PanOptic’s breadth of the field of view gives me a much greater visual appreciation of retinal structures. It increases my confidence during an exam.”

— LARRY CULPEPPER, MD,
BOSTON, MA



PanOptic™ Ophthalmoscope

The revolutionary PanOptic provides easier entry into the eye, together with a dramatically wider field of view and more magnification to help diagnose and manage diabetic and hypertensive retinopathy, papilledema, and retinal hemorrhage.

- Fast, easy entry into undilated pupils
- 5x larger view of the fundus vs. standard scopes in an undilated eye
- 25° field of view vs. the standard 5° field of view
- Greater working distance improves comfort for both practitioner and patient
- Detect corneal abrasions with cobalt blue filter and magnifying lens
- Compatible with all existing Welch Allyn 3.5v power sources
- Free diagnostic wall chart available

#11820

#11820 PanOptic™ Ophthalmoscope with Cobalt-Blue Filter and Add-On Corneal Viewing Lens (Head Only)

#11810 PanOptic™ Ophthalmoscope (Head Only)

3.5v Standard Coaxial Ophthalmoscope

- Halogen HPX™ lamp provides 30% whiter, brighter light for true tissue color and long-lasting performance
 - Coaxial optical system produces shadow-free spot viewing
 - Polarizing filter virtually eliminates corneal reflection
 - Sealed optics keep out dust and dirt
 - Detect corneal abrasions with cobalt blue filter (#11720 only)
- #11720** 3.5v Coaxial Ophthalmoscope (Head Only)
#11710 3.5v Standard Ophthalmoscope (Head Only)

#71910

3.5v Lithium Ion Rechargeable Handle

Twice the power, half the weight.

- Smaller diameter and half the weight of 3.5v Ni-Cad handles
- 120 minutes of operation between charges (compared to 60 minutes for standard Ni-Cad handles)
- Charges in new Universal Desk Charger or with removable AC adapter
- Compatible with all Welch Allyn 3.5v instrument heads
- Low battery indicator light

#71910 Handle alone, for use with Universal Desk Charger

#71900 Handle with Removable AC Adapter



3.5v Diagnostic Otoscopes

Featuring Halogen HPX™ illumination and fiber optics to the distal tip for easy examinations.

- Halogen HPX™ lamp provides 30% more light output for true tissue color and long-lasting performance
- Fiber optics project cool light to the distal tip with no reflections or obstructions, for clear viewing of the tympanic membrane
- Simply detach the head for high-intensity throat illumination or all-purpose penlight (models 20000 and 20001)
- Sliding viewing lens enables instrumentation under magnification
- Sealed system allows for pneumatic otoscopy and accurate diagnosis of otitis media

#20000 3.5v Diagnostic Otoloscope with Throat Illuminator (Head Only)

#25020 3.5v Diagnostic Otoloscope (Head Only)

#20001 3.5v Diagnostic Otoloscope with Throat Illuminator and Insufflator Bulb (Head Only)



2.5v PocketScope™ Set

Have all the accuracy of Welch Allyn diagnostics right in your lab coat pocket.

- 2.5v PocketScope Ophthalmoscope with consistent, long-lasting halogen illumination for true tissue color
- 2.5v Otoloscope with built-in throat illuminator gives you two instruments in one, with halogen illumination and fiber optics for cool light
- Adjustable rheostat to control light intensity
- Powered by AA batteries
- Available in either hard or soft case

#92820 PocketScope™ Set with Hard Case

#92821 PocketScope™ Set with Soft Case

Professional Stethoscopes

Excellent acoustics and superior design make this full-range instrument the first choice for professionals who need auscultation accuracy and all-day comfort.

- Meets all diagnostic needs, from detecting Korotkoff sounds of blood pressure to analyzing subtle heart murmurs
- Stainless steel chestpiece intensifies body sounds transmitted through tissues

#5079-135 Professional Adult, Double-head Stethoscope, 28" Tubing Shown in black. Available in a variety of colors.



#5079-135

Disposable Uniscope

The best protection against cross-contamination.

- Reduces risk of cross-contamination with the accuracy of conventional stethoscopes
- The safest choice for isolation, HIV and heavy fluid cases

#17462 Adult, Single-Head Disposable Stethoscope, Dark Blue

#17462P Pediatric, Single-Head Disposable Stethoscope, Light Blue




#17462

Lightweight Stethoscopes

Economical and comfortable for general use.

- Quality stethoscope designed for blood pressure measurement and general auscultation
- Lightweight and flexible single-channel tubing provides comfortable use over extended periods

#5079-73 Lightweight Double-Head Stethoscope Velvet black. Available in a variety of colors.



#5079-73

Durable and Disposable One-Piece Cuffs and Inflation Systems

Never remove, clean or replace a bladder.

- Cuffs are light, durable and latex-free
- Available in one- and two-tube designs for electronic NIBP monitors and sphygmomanometers
- Connectors attached to fit virtually every brand of NIBP monitors



See page 31 for complete listing of cuffs and order numbers.


Tycos® 767-Series Wall and Mobile Aneroids

Clinical-grade accuracy and mercury-free design make the 767-Series one of the best Tycos® aneroids ever.


- Certified accurate to ± 3 mmHg
- Mercury-free and latex-free for safety
- Ruby jeweled movement for long life
- Laser-engraved dial face for accuracy
- Large, easy-to-read dial face (6.25" / 15 cm)
- Convenient built-in cuff storage basket
- Adjustable-height, 360° swivel stand for clear viewing
- Lifetime warranty

#7670-01 Wall Aneroid with Durable One-Piece Cuff (Adult)

#7670-04 Mobile Aneroid with 4-Leg Mobile Stand, Durable One-Piece Cuff (Adult), and Base Weight



#7670-01



#5091-41

Tycos® 509-Series Wall and Mobile Aneroids


All the accuracy and reliability you expect from Tycos®—at an affordable price.

- Mercury-free and latex-free for safety
- Precision jeweled bearings resist wear for consistently accurate readings
- Thermally treated diaphragm maintains calibration
- Large, easy-to-read dial face
- Certified accurate to ± 3 mmHg
- Lifetime warranty

#5091-41 Mobile Aneroid with 4-Leg Stand and Durable One-Piece Cuff (Adult)

#5091-60 Mobile Aneroid with 5-Leg Stand and Durable One-Piece Cuff (Adult)

#5091-38 Wall Aneroid with Durable One-Piece Cuff (Adult)



#5098-27

Tycos® TR-1 Hand Aneroid

Innovative trigger-release design for easy operation.

- Certified accuracy to ± 3 mmHg
- Trigger deflation valve for smooth, one-finger operation
- Latex-free for safety
- Durable one-piece cuff
- Polyester zip case
- 10-year warranty

#5098-27 TR-1 Hand Aneroid with Durable One-Piece Cuff (Adult)

Welch Allyn® DuraShock™ Model DS45 Integrated Aneroid

The first gear-free, shock-resistant aneroid sphygmomanometer with integrated cuff.

- Meets AAMI shock-resistant standard for enhanced durability
- Laser-engraved dial face for accuracy
- Gauge rotates 360° for easy viewing and snaps directly into cuff for quick cuff change
- Lightweight and latex-free for safety
- Durable one-piece cuff
- 10-year warranty

#DS45-11 Sphygmomanometer Complete

#DS45-11C Sphygmomanometer with Case



#DS45-11



#5098-02

Tycos® Classic Hand Aneroid

Unsurpassed reliability in a traditional, handheld design.

- Certified accuracy to ± 3 mmHg
- Latex-free for safety
- Ruby jeweled movement for long life
- Laser-engraved dial face for accuracy
- Durable one-piece cuff
- Protective zipper case
- Lifetime warranty

#5098-02 Classic Hand Aneroid with Durable One-Piece Cuff (Adult)

Tycos® Classic Pocket Aneroid

All of the features of the Tycos Classic Hand Aneroid in a pocket design. The gauge attaches to the cuff, which prevents the cuff and gauge from being easily separated.

#5090-02 Classic Pocket Aneroid with Durable One-Piece Cuff (Adult)



Welch Allyn® Vital Signs Monitor 300 Series

For simple and accurate spot-check vital signs assessment or continuous monitoring, the newly redesigned Vital Signs Monitor 300 Series is the definition of “advanced monitoring made affordable.”

- Standard non-invasive blood pressure (NIBP), pulse rate and MAP
- Optional – SureTemp® Plus Thermometry*, motion tolerant Nellcor® Pulse Oximetry and integrated thermal printer
- Supports adult, pediatric and neonatal patients
- Simple to use – new LCD display facilitates menu navigation
- Includes automatic blood pressure mode, programmable alarms, built-in memory capability, and intuitive icons
- Lightweight and portable, with a variety of mounting options

- #53NTP-E1 Vital Signs Monitor complete with NIBP, Nellcor® Pulse Oximetry, SureTemp® Plus Thermometry and Printer
- #53000-E1 Vital Signs Monitor with NIBP (includes Pulse Rate and MAP)
- #530T0-E1 Vital Signs Monitor with NIBP and SureTemp Plus Thermometry
- #53N00-E1 Vital Signs Monitor with NIBP and Nellcor® Pulse Oximetry
- #53NT0-E1 Vital Signs Monitor with NIBP and Nellcor® Pulse Oximetry, and SureTemp Plus Thermometry
- #008-0835-00 Mobile Stand

*See page 14 for full details on Welch Allyn SureTemp® Plus Thermometry

Welch Allyn gives you automated vital signs cost-efficiently for procedures, minor surgeries, treatments, and routine assessments.

Spot Vital Signs®


Take vital signs in seconds with superior technology that spot checks easily, accurately and affordably.

- Automated measurement of non-invasive blood pressure (NIBP), pulse rate, temperature, and Pulse Oximetry in 30 seconds
- All configurations include: NIBP, Pulse Rate, and MAP
- Optional—SureTemp® Thermometry and motion-tolerant Nellcor® Pulse Oximetry
- Flexible configurations, yet fully upgradeable to meet your future needs
- Lightweight and portable with a variety of mounting options

- #42NTB-E1 Spot Vital Signs with NIBP, Nellcor® Pulse Oximetry, SureTemp® Thermometry
- #4200B-E1 Spot Vital Signs with NIBP (includes Pulse Rate and MAP)
- #420TB-E1 Spot Vital Signs with NIBP and SureTemp® Thermometry
- #42NOB-E1 Spot Vital Signs with NIBP, Nellcor® Pulse Oximetry
- #4200-60 Mobile Stand



NOW WITH MOTION TOLERANT SpO₂



#42NTB-E1



#42NTB-E1 & #4200-60



#53NTP-E1



“We love the convenience. Being able to automatically collect all the vital signs in a short amount of time really frees up the nurses.”

— LAURIE MCGRAW, RN,
MINONG, WI

To meet the needs of busy clinicians, Welch Allyn offers you all the modalities for temperature-taking, for all types of patients.



SureTemp® Plus 692 Thermometer

Welch Allyn's newly redesigned SureTemp Plus combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. You can take fast, four-second oral, 10-second pediatric axillary, 15-second adult axillary and 10-second rectal temperatures.

- New ergonomic design – comfortable to hold, easy to use, and very durable
- Interchangeable, removable probe wells to reduce the risk of cross-contamination
- Large backlit LCD display with location identification capability
- Includes a 60-second pulse timer for counting pulses and respirations
- Your choice of security alarm systems to help prevent theft
- Takes approximately 5000 readings on 3 AA batteries

#01692-200 SureTemp Plus 692


Thermometer, Oral Probe with Oral Probe Well, Wall Mount and Security System

#01692-700 SureTemp Plus 692 Thermometer, Oral Probe with Oral Probe Well, Wall Mount and Rolling Stand

#02892-000 Rectal Probe and Well Kit (4 ft)

#02893-000 Oral Probe and Well Kit (4 ft)

#05031-110 Thermometer Probe Covers (box of 10,000)



#01692-200



BRAUN ThermoScan. Pro3000

For patients and situations that require fast, accurate, non-invasive temperature measurement, the Braun Pro3000 infrared ear thermometer measures body temperature accurately in just one second.

- Sensor scans eight temperatures in one second, displaying the highest reading
- Built-in safeguards virtually eliminate risk of cross-contamination
- Soft, comfortable probe covers keep probe tip clean for accuracy
- Easy-to-read LCD in Fahrenheit or Celsius
- Takes approximately 5,000 readings on two lithium batteries

#03000-200 Braun ThermoScan® Pro 3000

#03000-300 Braun ThermoScan® Pro 3000 with Wall Mount

#03000-400 Braun ThermoScan® Pro 3000 with Security Roll Stand

#05074-005 Probe Covers (Box of 5000)



#03000-200

#05074-005

“After just one year, the increases in staff productivity made possible by Welch Allyn’s Flexible Monitoring® achieved reductions in patient transfers and I.o.s., resulting in significant savings for our hospital.”


– AILEEN EVANS, VICE PRESIDENT, PATIENT CARE SERVICES
METROPLEX HOSPITAL

LS-150 Halogen Exam Light

Totally re-engineered for performance and accuracy, our premium-quality exam light delivers brilliant color and convenience.

- Color-corrected 50-watt lamp delivers superior light quality for more accurate color rendition
- Specially designed light head keeps luminaire cool to the touch, even during multiple exams
- Reduced base profile frees up valuable exam room floor space
- Encased, flexible 24-inch arm enables one-hand positioning which is not possible with standard goose neck units
- Optional caster base for portability

#44500 LS-150 with Floor Stand
#44510 LS-150 with Wall Mount
#44350 Caster base accessory



#44350

LS-135 Halogen Exam Light

A unique blend of intense light output, mechanical stability, and affordability.

- Long-lasting, 35-watt halogen lamp produces significantly whiter, brighter light than ordinary goose neck units
- Dual-intensity control switch allows instant adjustment of light output
- Molded Lexan® head does not conduct heat, keeping luminaire cool to the touch at all times
- Encased, flexible 24-inch arm enables one-hand positioning which is not possible with standard goose neck units
- Optional caster base for portability

#44300 LS-135 with Floor Stand
#44310 LS-135 with Wall Mount
#44350 Caster base accessory



#76600

Professional PenLite

For quick and effective evaluation of the pupillary reflex.

- Bright, white halogen light for true tissue color rendition
- Durable solid brass construction

#76600 Professional PenLite



Welch Allyn Flexible Monitoring provides the freedom to monitor patients anywhere, anytime, lowering costs and increasing productivity.

Clinicians know the reality of hospital care. Keeping patients in high-acuity settings costs money. Moving the patient to the monitor—and back—takes time and hampers your efficiency.

Welch Allyn knows how you work. And Flexible Monitoring is one of the ways we're helping you work better. It makes any bed a monitored bed in seconds by letting you bring the monitor to the patient, and by giving you the freedom to continue monitoring that patient anywhere, anytime.

Flexible monitoring lets you care for patients efficiently and effectively, and in the process lowers the cost of monitoring and manages your time and resources better than ever.



Propaq® CS Vital Signs Monitor




Small, light, and rugged, it's the perfect instrument for bedside monitoring of vital signs.

- Large color display may be read across the room with color-differentiated parameters
- Custom selection of up to four viewable waveforms
- Adult, pediatric, and neonatal modes for customized monitoring of specific patient populations
- Optional—Wireless connectivity to the Acuity® and Acuity® LT Central Station allows truly flexible monitoring
- Three-year warranty

#01031633 Propaq CS, Model 242, Monitors 3-Lead ECG, Respiration, NIBP, Motion-Tolerant SpO₂, and 2 Channels of Temperature

#01031647 Propaq CS, Model 244, Monitors 3-Lead ECG, Respiration, NIBP, 1 Invasive BP, Motion-Tolerant SpO₂, Expansion Module with Printer, 2 Channels of Temperature, and Sidestream CO₂

#01031659 Propaq CS, Model 246, Monitors 3-Lead ECG, Respiration, NIBP, 2 Invasive BP, Motion-Tolerant SpO₂, Expansion Module with Printer, 2 Channels of Temperature, and Dualstream CO₂



Propaq Encore® Monitor

The gold standard for transport monitoring.

- Long battery life protects against data loss during transport and procedures
- Motion-tolerant SpO₂ and NIBP for accurate, consistent readings in challenging environments
- Low cost of ownership

#01020821 Propaq Encore, Model 202-EL, Monitors 3 Lead ECG, Respiration, NIBP, Motion-Tolerant SpO₂, and 2 Channels of Temperature


#01020834 Propaq Encore, Model 204-EL, Monitors 3 Lead ECG, Respiration, NIBP, 1 Invasive BP, Motion-Tolerant SpO₂, Expansion Module with Printer, 2 Channels of Temperature, and Sidestream CO₂

#01020846 Propaq Encore, Model 206-EL, Monitors 3 Lead ECG, Respiration, NIBP, 2 Invasive BP, Motion-Tolerant SpO₂, Expansion Module with Printer, 2 Channels of Temperature, and Dualstream CO₂






Acuity® Central Monitoring Station



Monitor more patients for less money.

- Monitor wireless and bedside devices on the same system
- Customized map of care unit shows each patient's status with easy-to-read, color-coded graphics
- Remote control of bedside monitors from the central station increases staff efficiency and enhances patient safety
- Easily expands to meet your facility's growing needs

Each Acuity and Acuity LT System is configured to your exact needs. Contact your Welch Allyn sales representative to plan and order the ideal configuration for your facility's size, patient load and floor layout.



Micropaq®

The convenience and efficiency of a patient monitor that's wireless and wearable. Ideal for monitoring vital signs of ambulating patients as well as administering medications at the bedside.

- Displays heart rate, one or two ECG leads, SpO₂ numerics and pulse bar
- Continues monitoring even when out of telemetry range
- Up to 25-hour battery life when fully charged
- Connects wirelessly to Welch Allyn's Acuity and Acuity LT Central Monitoring Stations
- Alarms may be suspended and silenced locally at the monitor, or remotely at the central station

- #01080004 ECG Monitoring Only
- #01080005 ECG and Masimo® Pulse Oximetry
- #01080006 ECG and Nellcor® Pulse Oximetry

#01080006



Acuity® LT Central Monitoring Station

Introducing the Acuity LT Central Monitoring Station, the affordable, self-contained, self-deployable, central monitoring solution. Used in conjunction with Welch Allyn's Propaq series and new Micropaq ambulatory monitors (see above), Acuity LT offers continuous monitoring for up to 12 patients — an excellent solution for providing surveillance in specific areas, rather than transferring patients to high-cost areas such as the ICU.

- Receives, displays and stores vital signs information for up to 12 patients from either a wireless or hardwired system of patient monitors
- Remote control of ambulatory or bedside patient monitors from the central station
- Twenty-four hour full disclosure and life-threatening arrhythmia analysis (Asystole, VFib and VTach alarms)
- Simplified installation allows your facility the cost-saving option of completing your own installation





#972212

Welch Allyn AED 20™

The unique, simple-to-operate AED 20 provides basic AED function with optional “Lead II” monitoring using traditional snap-based electrodes. Optional manual defibrillation (2-360 joules) is also available in this compact device. The versatile AED 20 is the first choice for health-care or EMS environments with multiple levels of first responders.

- Choice of rechargeable battery for frequent users or a long life non-rechargeable battery for infrequent rescuers
- Escalating energy (200j, 300j, 360j), per AHA guidelines. Energy protocol can be changed to meet specific needs
- Optional ECG monitoring and optional passcode protected manual defibrillation
- Data card transfer of critical patient data and direct connect printing to PCL5 printers provide efficient management of data


#972212 AED 20 with Rechargeable NiMH Battery and Charging Station

#972214 AED 20 with ECG Waveform and Single Channel Monitoring

#972216 AED 20 with ECG, Monitoring and Manual Defibrillation Mode

Welch Allyn AED 10™

Designed for easy operation and quick deployment by the infrequent user, the AED 10 uses the simple “1-2-3” approach – a proven rescue protocol used throughout the world. The AED 10 is the simple solution for quick and easy response by minimally trained personnel.



#970300

- Text and voice instructions guide the user through the rescue
- Color-coded pads and illustrations ensure proper pad placement
- Escalating energy (200j, 300j, 360j), per AHA guidelines. Energy protocol can be changed to meet specific needs
- High performance battery delivers five years of standby performance or over 90 discharges at 360 joules
- Internal memory stores event data including ECG waveforms, time/date stamps, shock delivery and system actions
- Compatible with major manufacturer defibrillator pads

970300 AED 10


970302 AED 10 with Soft Shell Carrying Case

Welch Allyn PIC 50™

The PIC 50 is a comprehensive Portable Intensive Care Unit, with up to 12 patient treatment parameters and sophisticated communication options. Designed for high acuity settings, the PIC 50 boasts multi-parameter monitoring and defibrillation capabilities in one compact, lightweight and durable package.

- 3 or 5 lead ECG monitoring, multi-parameter defibrillator, recorder and SpO₂
- Options include: Full diagnostic 12 lead ECG with interpretation, AED, ETCO₂, non-invasive blood pressure, two invasive blood pressure lines, external pacing, and temperature, all in a single 10 lb. box
- Large active matrix color display with color-coded parameters that can be viewed from virtually any angle
- Includes high performance, lightweight NiMH battery to deliver more than four hours of continuous performance from a single charge
- A breadth of power options and battery support systems are also available to meet virtually any need
- Optional data card transfer, fax, GSM cellular and other sophisticated communications options are also available

PART NUMBERS: Contact your local Welch Allyn Hospital sales representative to plan and order the ideal PIC 50 configuration for your facility, patient and treatment needs.



#973039C00

Welch Allyn PIC 30™

Designed for the crash-cart and other manual defibrillation environments, the PIC 30 monitor/defibrillator is the economical choice for your patient care department.

- Radiant color display provides excellent visibility
- Annotating recorder captures all vital data on 50mm printer
- Full manual defibrillation (2 joules to 360 joules) and sync capabilities
- Optional external pacing and shock advisory
- Also available with monochrome display

#973039C00 PIC 30 with Color Screen

#973039COA PIC 30 with Shock Advisory Feature

#973039CP0 PIC 30 with External Pacing

#973039CPA PIC 30 with External Pacing and Shock Advisory

Welch Allyn PIC 40™

A multi-parameter monitor, AED, and defibrillator with external pacing, the PIC 40 is the perfect solution for procedure areas.

- Radiant color display provides waveforms, patient data and defibrillation functions that can be viewed from virtually any angle
- 3 or 5 lead ECG and SpO₂
- Optional parameters include: NIBP and external pacing
- Manual defibrillator provides full range of energy settings (2 joules to 360 joules)
- AED is configurable and follows AHA escalating energy guidelines (200j, 300j, 360j)

#973092 PIC 40

#973093 PIC 40 with NIBP

#973094 PIC 40 with NIBP and External Pacing

#973095 PIC 40 with External Pacing



#973092

Atlas™ Monitor

A cost-effective monitor with all the parameters for conscious sedation procedures.



- Lean, practical, hard-working monitor with configuration options to meet your individual needs
- All monitors include ECG, motion-tolerant Nellcor® Pulse Oximetry and NIBP. Optional parameters include temperature, respiration, ETCO₂ and integrated printer.
- Bright, crisp waveforms and multi-colored LED displays for across-the-room visibility
- Simple, straightforward operation
- Two-year warranty

#623NP-E1 Atlas™ with ECG, Nellcor® Pulse Oximetry, NIBP, ETCO₂, Temperature, Respiration, Battery, RS232, and Printer


#622N0-E1 Atlas™ with ECG, Nellcor® Pulse Oximetry, NIBP, Temperature, Respiration, Battery, and RS232

#622NP-E1 Atlas™ with ECG, Nellcor® Pulse Oximetry, NIBP, Temperature, Respiration, Battery, RS232, Printer

#6200-30 Mobile Stand



#623NP-E1




SPECIALTY PROCEDURES

“Welch Allyn gives me a level of confidence that no one else comes close to. No matter which one of their instruments I'm using, I never worry about reliability, accuracy or dependability. I know I can rely on Welch Allyn devices to work consistently and work well, time after time.”

— KATHLEEN KENDRICK, RN, MBA, CCRN, NURSE MANAGER
COMMUNITY GENERAL HOSPITAL, SYRACUSE, NY

“I have enjoyed using Welch Allyn’s disposable specula for over five years. The visualization is far superior to the old metal specula with goose-neck lamp. I had the opportunity to use the newest speculum with the attached sheath... it was a pleasure. Thank you Welch Allyn.”

— LINDA S. MCCARTNEY, MD, FAAEM, CHICAGO, IL



KleenSpec® Vaginal Specula and Illumination System

Welch Allyn’s new KleenSpec® with attached sheath reduces the risk of illuminator contamination.

- Preattached sheath easily pulls down over the Welch Allyn illumination system
- The deployed sheath protects the light cord from cross-contamination
- Available in small, medium and large-size specula
- Unique Welch Allyn high-quality molded acrylic for patient comfort and reliability
- Specula also available without sheath

#78810 KleenSpec® Illumination System with On/Off Switch and Transformer



- #58000S KleenSpec® Small Disposable Specula with Attached Sheath (Case of 80)
- #58001S KleenSpec® Medium Disposable Specula with Attached Sheath (Case of 80)
- #58004S KleenSpec® Large Disposable Specula with Attached Sheath (Case of 60)
- #58000 KleenSpec® Small Disposable Specula (Case of 100)
- #58001 KleenSpec® Medium Disposable Specula (Case of 100)
- #58004 KleenSpec® Large Disposable Specula (Case of 80)



AVAILABLE IN 2004

KleenSpec® Illumination Storage System

Store your KleenSpec® illumination system in this convenient storage unit that directly mounts wherever needed. Organizes the illuminator, cords, and extra lamps.

#78820 KleenSpec® Illumination Storage System



Fiber Optic Laryngoscope Handles

Bright, white HPX™ illumination offers a clear, color-corrected view.

- Patented HPX™ technology blends halogen and xenon gases under pressure to produce the brightest, whitest light available, for true tissue color and consistent, long-lasting illumination
- Single lamp cartridge can be used across all sizes and handles
- Up to 40% lighter than our previous design, for increased responsiveness and maneuverability
- Knurled finish ensures secure grip

Lightweight Battery-Powered Fiber Optic Handles

Model	Description
#60813	Medium, uses two C-size batteries
#60814	Penlight, uses two AA-size batteries
#60815	Stubby, uses two AA-size batteries

Fiber Optic Laryngoscope Blades

- Removable, repolishable light pipes and free repolishing service help blades last longer
- Damaged light pathways can be repaired at a fraction of the cost of replacement blade
- Single-piece stainless steel construction ensures blade integrity

MacIntosh Fiber Optic Blades

Model	Blade Size	Overall Length	Inside Base to Tip
#69061	#1	89 mm	63 mm
#69062	#2	108 mm	82 mm
#69063	#3	126 mm	101 mm
#69064	#4	160 mm	135 mm

English MacIntosh Fiber Optic Blades

Model	Blade Size	Overall Length	Inside Base to Tip
#69211	#1	95 mm	70 mm
#69212	#2	116 mm	90 mm
#69213	#3	135 mm	110 mm
#69214	#4	156 mm	130 mm

Miller Fiber Optic Blades

Model	Blade Size	Overall Length	Inside Base to Tip
#68065	#00	8 mm	37 mm
#68060	#0	8 mm	53 mm
#68061	#1	8 mm	79 mm
#68062	#2	10 mm	132 mm
#68063	#3	11 mm	172 mm
#68064	#4	13 mm	182 mm

LS-200™ Procedure Light

Bring light to all your procedures with high-quality, low-cost illumination.

- Superior maneuverability covers a wide range of positions without drifting or repositioning
- 24" working distance with 2,500 foot-candles light output and 3,200 K color temperature enhances visualization of the operative site and true tissue color
- Three-lamp design ensures continuous, shadow-free illumination
- 7" spot size with feathered edge definition and 12" depth of field reduces eye strain

#44200 LS-200 Procedure Light

#44250 Mobile Stand

#44210 Wall Mount Bracket

Solid State Portable Headlight

Provides illumination precisely where it is needed, with cool, bright light that ensures comfort for both clinician and patient.

- Weighs just 4 ounces — one of the most comfortable and lightweight headlamps available
- Fits inside a shirt pocket and allows the user complete freedom



of mover(2. ide a shirt pock907001 Tc0.027.f mover(2. ide a Tj/F3 1 Tf2.9072 0 Tm-pliW)7cool,) isehf t, res continuoT 0 -ht tT*(operative site ad reduces eye strain)Tj/F1 1 Tf46.9072 -1.444607D-0.0002 Tc03.24 144210Prt PowpocSTwce of tChargeite ad reduces 3.24 l4and

Consumables Reference Chart

Product Name	Part No.	Ref. Page	Product Affiliation
Eartips	#52434	8	524 Series KleenSpec® Disposable Otoscope Specula, 4 mm (1,000/bag, 10 bags/case)
	#52432	8	524 Series KleenSpec® Disposable Otoscope Specula, 2.5 mm (1,000/bag, 10 bags/case)
	#52401	6	524 Series Disposable Otoscope Specula Dispenser with Storage Compartment (case of 10)
Thermometry Probe Covers	#05031-110	14	SureTemp® Disposable Probe Covers (10,000 covers)
	#05031-105	14	SureTemp® Disposable Probe Covers (5,000 covers)
	#05074-005	15	Braun Disposable Probe Covers for Pro3000 (case of 5,000)
Paper Supplies	#7052-25	13	Vital Signs Monitor Printer Paper, 1 package (6 rolls)
	#6200-40	24	Atlas™ Monitor Printer Paper, 1 case (25 rolls)
Batteries	#72200	5, 6	Replacement Ni-Cad Rechargeable Battery (black) for Desk/Well Charger Handle
	#71960	7, 8	Replacement Lithium-Ion Battery
Lamps	#03000	7	Ophthalmoscope #11710
	#04900	7	Ophthalmoscopes #11720, 11730, 11735
	#03800	7	PanOptic™ Ophthalmoscopes #11810, 11820
	#03100	8	Otosopes #20000, 25020, 21700, 20200, 3.5v Handle Adapter #73500
	#04430	16	35w Halogen Lamp for LS-135 Exam Light #44300 and 44310
	#04450	16	50w Lamp for LS-150 Exam Light #44500 and 44510
	#08800-U	26	4.6v Halogen Lamp for Vaginal Specula Illumination System #78810
	#06000	27	Halogen HPX™ Lamp for Fiber Optic Laryngoscope Handles
Electrodes	#04200	28	Halogen Lamp for Exam Light II™ #48600 and Exam Light III™ #48700
	#06400	28	20w Lamp for LS-200™ Procedure Light #44200
	#6200-11	25	ECG Electrodes for Atlas Monitor (pack of 30)
Vaginal Specula	#58000	26	KleenSpec® 580 Series Disposable Vaginal Specula, Small (case of 100)
	#58001	26	KleenSpec® 580 Series Disposable Vaginal Specula, Medium (case of 100)
	#58004	26	KleenSpec® 580 Series Disposable Vaginal Specula, Large (case of 80)
	#58000S	26	KleenSpec® Small Disposable Specula with Attached Sheath (case of 80)
	#58001S	26	KleenSpec® Medium Disposable Specula with Attached Sheath (case of 80)
	#58004S	26	KleenSpec® Large Disposable Specula with Attached Sheath (case of 60)
Miscellaneous	#76700	5, 6	3.5v Locking Collars for Wall Transformers
	#11870	7	PanOptic™ Patient Eye Cups (package of 5)
	#001855	22, 23	Defibrillator Pads, Adult (box of 10)
	#52640	28	Disposable Sheaths for Exam Light II™ and Exam Light III™ (#48600/48700)

Durable One-Piece Cuffs and Inflation Systems

MONITOR STYLE, ONE TUBE, SCREW CONNECTOR

Size	Circum.	Single	5 Pack*
Newborn	7.7-10.5 cm	5082-201-3	5082-81-3
Infant	9.8-13.3 cm	5082-202-3	5082-82-3
Small Child	12.4-16.8 cm	5082-203-3	5082-83-3
Child	15.8-21.3 cm	5082-204-3	5082-84-3
Child Print	15.8-21.3 cm	5082-204P-3	5082-84P-3
Small Adult	20.0-27.0 cm	5082-205-3	5082-85-3
Adult	25.3-34.3 cm	5082-206-3	5082-86-3
Adult Long	25.3-34.3 cm	5082-206L-3	5082-86L-3
Large Adult	32.1-43.4 cm	5082-207-3	5082-87-3
Large Adult Long	32.1-43.4 cm	5082-207L-3	5082-87L-3
Thigh	40.7-55.0 cm	5082-208-3	5082-88-3

Special Notes: To order multi-packs with one of each size cuff (excluding Child Print and Long Cuffs), order 5082-241-3.
 *To order cuffs with bayonet connectors (#5082-184), add "HP" suffix to part number (e.g., 5082-81-3HP).
 *To order cuffs with locking connectors (#5082-182), add "MQ" suffix to part number (e.g., 5082-81-3MQ).

MONITOR STYLE, TWO TUBES, SCREW CONNECTORS

Size	Circum.	Single	5 Pack*
Infant	9.8-13.3 cm	5082-202-4	5082-82-4
Small Child	12.4-16.8 cm	5082-203-4	5082-83-4
Child	15.8-21.3 cm	5082-204-4	5082-84-4
Child Print	15.8-21.3 cm	5082-204P-4	5082-84P-4
Small Adult	20.0-27.0 cm	5082-205-4	5082-85-4
Adult	25.3-34.3 cm	5082-206-4	5082-86-4
Adult Long	25.3-34.3 cm	5082-206L-4	5082-86L-4
Large Adult	32.1-43.4 cm	5082-207-4	5082-87-4
Large Adult Long	32.1-43.4 cm	5082-207L-4	5082-87L-4
Thigh	40.7-55.0 cm	5082-208-4	5082-88-4

Special Notes: To order multi-packs with one of each size cuff (excluding Child Print and Long Cuffs), order 5082-241-4.
 *To order cuffs with bayonet connectors (#5082-184), add "HP" suffix to part number (e.g., 5082-82-4HP).
 *To order cuffs with locking connectors (#5082-182), add "MQ" suffix to part number (e.g., 5082-82-4MQ).

SPHYG STYLE, ONE TUBE, LUER LOCK CONNECTOR

Size	Circum.	Single	5 Pack
Newborn	7.7-10.5 cm	5082-201-1	5082-81-1
Infant	9.8-13.3 cm	5082-202-1	5082-82-1
Small Child	12.4-16.8 cm	5082-203-1	5082-83-1
Child	15.8-21.3 cm	5082-204-1	5082-84-1
Child Print	15.8-21.3 cm	5082-204P-1	5082-84P-1
Small Adult	20.0-27.0 cm	5082-205-1	5082-85-1
Adult	25.3-34.3 cm	5082-206-1	5082-86-1
Adult Long	25.3-34.3 cm	5082-206L-1	5082-86L-1
Large Adult	32.1-43.4 cm	5082-207-1	5082-87-1
Large Adult Long	32.1-43.4 cm	5082-207L-1	5082-87L-1
Thigh	40.7-55.0 cm	5082-208-1	5082-88-1

Special Note: To order multi-packs with one of each size cuff (excluding Child Print and Long Cuffs), order 5082-241-1.

SPHYG STYLE, TWO TUBES, LUER SLIP CONNECTORS

Size	Circum.	Single	5 Pack
Infant	9.8-13.3 cm	5082-202-2	5082-82-2
Small Child	12.4-16.8 cm	5082-203-2	5082-83-2
Child	15.8-21.3 cm	5082-204-2	5082-84-2
Child Print	15.8-21.3 cm	5082-204P-2	5082-84P-2
Small Adult	20.0-27.0 cm	5082-205-2	5082-85-2
Adult	25.3-34.3 cm	5082-206-2	5082-86-2
Adult Long	25.3-34.3 cm	5082-206L-2	5082-86L-2
Large Adult	32.1-43.4 cm	5082-207-2	5082-87-2
Large Adult Long	32.1-43.4 cm	5082-207L-2	5082-87L-2
Thigh	40.7-55.0 cm	5082-208-2	5082-88-2

Special Note: To order multi-packs with one of each size cuff (excluding Child Print and Long Cuffs), order 5082-241-2.

INFLATION SYSTEMS, TWO TUBES WITH BULB AND VALVE

Size	Circum.	Single
Infant	9.8-13.3 cm	5082-222-1
Small Child	12.4-16.8 cm	5082-223-1
Child	15.8-21.3 cm	5082-224-1
Child Print	15.8-21.3 cm	5082-224P-1
Small Adult	20.0-27.0 cm	5082-225-1
Adult	25.3-34.3 cm	5082-226-1
Adult Long	25.3-34.3 cm	5082-226L-1
Large Adult	32.1-43.4 cm	5082-227-1
Large Adult Long	32.1-43.4 cm	5082-227L-1
Thigh	40.7-55.0 cm	5082-228-1

CONNECTORS AND ADAPTERS FOR DURABLE AND DISPOSABLE CUFFS

5082-161	Metal screw-type adapter with barbed end for Welch Allyn and Critikon monitors
5082-164	Plastic male screw-type connector with female luer slip for Welch Allyn and Critikon monitor cuffs
5082-165	Plastic male locking luer slip connector with barbed end for Datascope monitors
5082-168	Plastic female locking luer slip connector with barbed end for sphygmomanometers and Datascope monitor cuffs
5082-169	Metal female luer slip adapter with barbed end
5082-174	Plastic screw-type adapter with barbed end for Welch Allyn and Critikon monitors
5082-176	Plastic male luer slip connector with barbed end for neonatal cuffs
5082-178	Plastic male locking adapter for Welch Allyn Protocol and GE Marquette monitors
5082-182	Plastic female locking connector for Welch Allyn Protocol and GE Marquette monitor cuffs
5082-184	Plastic male bayonet connector for HP/Siemens, Spacelabs, and Datascope monitors
5082-185	Metal female bayonet adapter for HP/Siemens, Spacelabs, and Datascope monitor cuffs

Disposable One-Piece Cuffs and Inflation Systems

MONITOR STYLE, ONE TUBE, SCREW CONNECTOR

Size	Circum.	Soft 5 Pack*	Vinyl 5 Pack*
Newborn	7.7-10.5 cm	5082-91-3	5082-91-8
Infant	9.8-13.3 cm	5082-92-3	5082-92-8
Small Child	12.4-16.8 cm	5082-93-3	5082-93-8
Child	15.8-21.3 cm	5082-94-3	5082-94-8
Small Adult	20.0-27.0 cm	5082-95-3	5082-95-8
Adult	25.3-34.3 cm	5082-96-3	5082-96-8
Adult Long	25.3-34.3 cm	5082-96L-3	5082-96L-8
Large Adult	32.1-43.4 cm	5082-97-3	5082-97-8
Large Adult Long	32.1-43.4 cm	5082-97L-3	5082-97L-8
Thigh	40.7-55.0 cm	5082-98-3	5082-98-8

Special Notes: To order multi-packs with one of each size cuff (excluding Long Cuffs), order 5082-241-7 for soft or 5082-241-19 for vinyl.
 *To order cuffs with bayonet connectors (#5082-184), add "HP" suffix to part number (e.g., 5082-91-3HP).
 *To order cuffs with locking connectors (#5082-182), add "MQ" suffix to part number (e.g., 5082-91-3MQ).

MONITOR STYLE, TWO TUBES, SCREW CONNECTORS

Size	Circum.	Soft 5 Pack*	Vinyl 5 Pack*
Infant	9.8-13.3 cm	5082-92-4	5082-92-9
Small Child	12.4-16.8 cm	5082-93-4	5082-93-9
Child	15.8-21.3 cm	5082-94-4	5082-94-9
Small Adult	20.0-27.0 cm	5082-95-4	5082-95-9
Adult	25.3-34.3 cm	5082-96-4	5082-96-9
Adult Long	25.3-34.3 cm	5082-96L-4	5082-96L-9
Large Adult	32.1-43.4 cm	5082-97-4	5082-97-9
Large Adult Long	32.1-43.4 cm	5082-97L-4	5082-97L-9
Thigh	40.7-55.0 cm	5082-98-4	5082-98-9

Special Notes: To order multi-packs with one of each size cuff (excluding Long Cuffs), order 5082-241-8 for soft or 5082-241-20 for vinyl.
 *To order cuffs with bayonet connectors (#5082-184), add "HP" suffix to part number (e.g., 5082-92-4HP).
 *To order cuffs with locking connectors (#5082-182), add "MQ" suffix to part number (e.g., 5082-92-4MQ).

MONITOR STYLE, ONE TUBE, LUER SLIP CONNECTOR

Size	Circum.	Soft 10 Pack
Neonate #1	3.3-5.6 cm	5082-101-1
Neonate #2	4.2-7.1 cm	5082-102-1
Neonate #3	5.4-9.1 cm	5082-103-1
Neonate #4	6.9-11.7 cm	5082-104-1
Neonate #5	8.9-15.0 cm	5082-105-1

Special Note: To order multi-packs with one of each size cuff, order 5082-241-9.

MONITOR STYLE, TWO TUBES, LUER SLIP CONNECTORS

Size	Circum.	Soft 10 Pack
Neonate #1	3.3-5.6 cm	5082-101-2
Neonate #2	4.2-7.1 cm	5082-102-2
Neonate #3	5.4-9.1 cm	5082-103-2
Neonate #4	6.9-11.7 cm	5082-104-2
Neonate #5	8.9-15.0 cm	5082-105-2

Special Note: To order multi-packs with one of each size cuff, order 5082-241-10.

SPHYG STYLE, ONE TUBE, LUER SLIP CONNECTOR

Size	Circum.	Soft 5 Pack	Vinyl 5 Pack
Newborn	7.7-10.5 cm	5082-91-1	5082-91-6
Infant	9.8-13.3 cm	5082-92-1	5082-92-6
Small Child	12.4-16.8 cm	5082-93-1	5082-93-6
Child	15.8-21.3 cm	5082-94-1	5082-94-6
Small Adult	20.0-27.0 cm	5082-95-1	5082-95-6
Adult	25.3-34.3 cm	5082-96-1	5082-96-6
Adult Long	25.3-34.3 cm	5082-96L-1	5082-96L-6
Large Adult	32.1-43.4 cm	5082-97-1	5082-97-6
Large Adult Long	32.1-43.4 cm	5082-97L-1	5082-97L-6
Thigh	40.7-55.0 cm	5082-98-1	5082-98-6

Special Note: To order multi-packs with one of each size cuff (excluding Long Cuffs), order 5082-241-5 for soft or 5082-241-17 for vinyl.

SPHYG STYLE, TWO TUBES, LUER SLIP CONNECTORS

Size	Circum.	Soft 5 Pack	Vinyl 5 Pack
Infant	9.8-13.3 cm	5082-92-2	5082-92-7
Small Child	12.4-16.8 cm	5082-93-2	5082-93-7
Child	15.8-21.3 cm	5082-94-2	5082-94-7
Small Adult	20.0-27.0 cm	5082-95-2	5082-95-7
Adult	25.3-34.3 cm	5082-96-2	5082-96-7
Adult Long	25.3-34.3 cm	5082-96L-2	5082-96L-7
Large Adult	32.1-43.4 cm	5082-97-2	5082-97-7
Large Adult Long	32.1-43.4 cm	5082-97L-2	5082-97L-7
Thigh	40.7-55.0 cm	5082-98-2	5082-98-7

Special Note: To order multi-packs with one of each size cuff (excluding Long Cuffs), order 5082-241-6 for soft or 5082-241-18 for vinyl.

INFLATION SYSTEMS, TWO TUBES WITH BULB AND VALVE

Size	Circum.	Soft Single	Vinyl Single
Infant	9.8-13.3 cm	5082-232-1	5082-232-6
Small Child	12.4-16.8 cm	5082-233-1	5082-233-6
Child	15.8-21.3 cm	5082-234-1	5082-234-6
Small Adult	20.0-27.0 cm	5082-235-1	5082-235-6
Adult	25.3-34.3 cm	5082-236-1	5082-236-6
Adult Long	25.3-34.3 cm	5082-236L-1	5082-236L-6
Large Adult	32.1-43.4 cm	5082-237-1	5082-237-6
Large Adult Long	32.1-43.4 cm	5082-237L-1	5082-237L-6
Thigh	40.7-55.0 cm	5082-238-1	5082-238-6