

SPECIALTY ITEMS

	Depth Adjustment	Height Adjustment
6088	24"—36"	26"—39"
	Width	Length
Armrests	6"	24"

Conveniently stores against wall out of the way of foot traffic.

Wheelchair Blood Drawing Station

- Wall-mounted for secure dependable performance
- Adjustable height and depth
- Armrest locks in place when in use, then flips out of the way for easy exit
- Comfortable 2" padding
- Easy-clean, antimicrobial, heavy-weight, vinyl upholstery
- 2 knobs lock chrome armrest into position
- Plastic work tray (optional)

	Width	Length	Height
6077	20 1/2"	39"	39 1/2"
	Width	Length	
Top	20"	32"	

Infant Blood Drawing Station

- All-welded stainless steel frame
- Top has 2 rims
- 2 safety straps
- Plastic work tray
- Removable middle padded section