

ELAN SERIES STOOLS

CONTEMPORARY, PREMIUM SEATING

FEATURES

- Color coordinated Natural Sand bases, column covers and backrest shrouds
- Premium, die-cast, powder-coated, aluminum base
- 5-leg design for greater stability
- Fast, pneumatic height adjustment
- Large, easy to locate "D-shaped" height lever
- Pneumatic lever is field changeable from left hand to right
- "Soft Feel" polyfoam padding
- Hooded, dual wheel, nylon casters
- All seats have solid plywood base
- 2 styles of contour backrest
- ABS shrouded backrest for easy cleaning
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- 23 standard color choices



Easy to clean backrest



"Soft Feel" foam provides two layers of comfort



	Seat Diameter	Height Range
2175-041	16"	17 ¹ / ₂ "–23"

5-Leg Pneumatic Stool with Lumbar-Shaped Backrest

- Dual wheel casters
- 23" powder-coated aluminum base
- Large chrome pneumatic lever
- Comfortable lumbar-shaped backrest with depth adjustment

	Seat Diameter	Height Range
2175	16"	17 ¹ / ₂ "–23"

5-Leg Pneumatic Stool

- Dual wheel casters
- 23" powder-coated aluminum base
- Large chrome pneumatic lever

	Seat Diameter	Height Range
2175-031	16"	17 ¹ / ₂ "–23"

5-Leg Pneumatic Stool with D-Shaped Backrest

- Dual wheel casters
- 23" powder-coated aluminum base
- Large chrome pneumatic lever
- Comfortable D-shaped backrest with depth adjustment

	Seat Diameter	Height Range
2175W-031	18"	16 ¹ / ₄ "–21 ³ / ₄ "

5-Leg Pneumatic Stool with Contour Seat and D-Shaped Backrest

- Dual wheel casters
- 23" powder-coated aluminum base
- Large chrome pneumatic lever
- Depth adjustable backrest
- Comfortable, contour seat & backrest
- Seamless upholstered seat & backrest

	Seat Diameter	Height Range
2175W-041	18"	16 ¹ / ₄ "–21 ³ / ₄ "

5-Leg Pneumatic Stool with Contour Seat and Lumbar-Shaped Backrest

- Same as 2175W-031 with lumbar-shaped backrest

Lumbar-shaped backrest for added support