

# ETA TREATMENT CABINETS

Clinton ETA (Easy-To-Assemble) Treatment Cabinets are designed with the needs of the professional in mind. Six models feature modules that ship easy and go together fast using only a Phillips screwdriver. Mix or match up to 12 colors of laminate for a traditional or unique look.

## SHIPPING ADVANTAGES

- Save money on shipping. Clinton ETA cabinets ship via UPS® or FedEx® Ground
- Single person delivery and assembly
- No inside delivery charges

### NOTE:


Large orders that are less expensive to be shipped by common carrier will be shipped assembled.

## FEATURES

- Pre-assembled modules
- All assembly holes pre-drilled and pre-threaded
- Easy-clean all laminate construction
- Tough chip-resistant PVC edges
- Modules are assembled with steel fasteners, not wood dowels
- 3-way, adjustable 110° door hinges
- Euro-style drawers with metal sides that hold up to 75 lbs.
- Recessed toe kick or dual wheel caster bases
- Height-adjustable drawer fronts
- Mix color choices for top sides and fronts
- Light color interior laminates
- Any damaged components are field changeable
- Assemble with only a Phillips screwdriver

### ORDERING INFORMATION:

When ordering cabinets you need to select top, front and side colors. If only one color is selected, that color will be used for the entire cabinet.



Clinton uses pre-drilled holes and pre-threaded fasteners that allow Clinton ETA cabinets to be assembled faster and easier.



	Length	Height	Depth
<b>8821</b>	24"	35"	18"

**Treatment Cabinet with 2 Doors & 1 Drawer**

- 1 adjustable shelf
- Euro-style drawer with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
<b>8822</b>	24"	35"	18"

**Treatment Cabinet with 2 Doors & 2 Drawers**

- 1 adjustable shelf
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
<b>8805</b>	24"	35"	18"

**Treatment Cabinet with 5 Drawers**

- Smooth glide tracks for easy open and close
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
<b>8922</b>	24"	34"	18"

**Mobile Treatment Cabinet with 2 Doors and 2 Drawers**

- Mounted on dual wheel swivel casters
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
<b>8921</b>	24"	34"	18"

**Mobile Treatment Cabinet with 2 Doors and 1 Drawer**

- Mounted on dual wheel swivel casters
- Euro-style drawer with metal sides that hold up to 75 lbs.
- All laminate construction

	Length	Height	Depth
<b>8950</b>	24"	34"	18"

**Mobile Treatment Cabinet with 5 Drawers**

- Mounted on dual wheel swivel casters
- Euro-style drawers with metal sides that hold up to 75 lbs.
- All laminate construction