

ALPHA SERIES MAT PLATFORMS

FEATURES

- All-welded, heavy duty steel frame finished in Gray
- Heavy-duty 1 1/2" x 3" steel legs
- Plywood top
- Solid six-leg support
- 18" fixed height—standard
- 1,000 lbs. weight capacity

	Width	Length	Height
230-47	4'	7'	18"
230-57	5'	7'	18"
230-68	6'	8'	18"

Upholstered Top Mat Platform

- 2" of firm foam padding covered with heavy weight vinyl with Perma-Block top coating
- Custom edge bumper guard helps to prevent top from tearing

	Width	Length	Height
231-47	4'	7'	18"
231-57	5'	7'	18"
231-68	6'	8'	18"

Upholstered Top Mat Platform with Adjustable Backrest

- 2" of firm foam padding covered with heavy weight vinyl with Perma-Block top coating
- Custom edge bumper guard helps to prevent top from tearing
- 6-position, all-steel, adjustable backrest
- 6' x 8' model has split 30" high backrest

	Width	Length	Height
232-47	4'	7'	18"
232-57	5'	7'	18"
232-68	6'	8'	18"

Mat Platform with Removable Top

- 2" mat attaches with Velcro® strips to the rimless subtop
- Mat is included