

HI-LO MAT PLATFORMS

	Width	Length	Height Range
251-47	4'	7'	23"-35"
251-57	5'	7'	23"-35"
251-68	6'	8'	23"-35"

Hand-Crank Hydraulic Hi-Lo Mat Platform

- All-welded steel frame with laminated finish sides
- Stable 4-leg lift system with levelers
- 2" thick firm padding
- Custom bumper strip protects top's upholstery
- Easy to operate hand-crank height adjustment
- Hand crank easily folds for storage
- 750 lbs. lift capacity

	Width	Length	Height Range
253-47	4'	7'	22"-34"
253-57	5'	7'	22"-34"
253-68	6'	8'	22"-34"

Hi-Lo Mat Platform

- All-welded steel frame
- Hands-free, foot peddle control
- Dependable single motor drive
- 110 volts, 60 hz., 2.7 amps. with overload protection
- Fully enclosed drive unit for safety
- Custom bumper strip protects top's upholstery
- 600 lbs. lift capacity

	Width	Length	Height Range
252-47	4'	7'	22"-34"
252-57	5'	7'	22"-34"
252-68	6'	8'	22"-34"

Hi-Lo Mat Platform with Removable Mat

- All-welded steel frame
- Hands-free, foot peddle control
- Dependable single motor drive
- 110 volts, 60 hz., 2.7 amps. with overload protection
- Fully enclosed drive unit for safety
- 600 lbs. lift capacity
- Removable 2" mat held in place with Velcro® strips
- Mat included