

PARALLEL BARS


FEATURES

- One piece stainless steel handrails with end bumpers
- Uprights are ultra-heavy, square, telescoping steel tubing
- All models, except wall mounted folding unit, feature width adjustment
- Platform models feature 1 1/2" thick, hardwood, plywood bases with ultra-durable, UV cured, satin finish
- Safety treads at both ends of platform models make it safer for patient to turn and exit wheelchair
- All models are for Physical Therapy use only

	Length	Width Range	Height Range
3-2007	7'	18"-28"	26"-39"
3-2010	10'	18"-28"	26"-39"
3-2012	12'	18"-28"	26"-39"

Platform Mounted Parallel Bars

- Stainless steel adjustable handrails
- 400 lbs. load capacity



PARALLEL BARS

	Length	Width	Height Range
3-3307	7'	28"	28"-42"

Wall Mounted Folding Parallel Bars

- Stainless steel handrails
- Requires only 14' of floor space when stored
- Folds up against the wall when not in use
- 400 lbs. load capacity

	Length	Width Range	Height Range
3-4007	7'	adjusts 6"/side	26"-39"
3-4010	10'	adjusts 6"/side	26"-39"
3-4012	12'	adjusts 6"/side	26"-39"

Floor Mounted Parallel Bars

- Stainless steel adjustable handrails
- Same features as platform-mounted model, except it is floor-mounted
- 400 lbs. load capacity



	Length	Width Range	Height Range
3-2106	10'	20"-34"	26"-39"

Bariatric Parallel Bars

- Height and width adjustment
- 3rd upright with through pin lock for extra support
- Extra width for bariatric patients
- 48" platform width
- 500 lbs. load capacity

	Length	Width Range	Height Range
3-3317	7'	16"-24"	22"-35"
3-3310	10'	16"-24"	22"-35"

Folding Parallel Bars

- Stainless steel handrails
- Unique double stabilizer bars stiffen the rails
- Spring-loaded plungers and sure lock positioning knobs
- 350 lbs. load capacity

PARALLEL BARS



	Length	Width	Height Range
3-5007	7'	20"	26"–39"
3-5010	10'	20"	26"–39"
3-5012	12'	20"	26"–39"

Hemiplegic Parallel Bars

- Floor mounted bars allow patients an uninterrupted bar to grasp with one hand
- 350 lbs. load capacity

-045

Optional Center Obduction Board

- Available with platform style only