


# CLASSIC SERIES STAIRCASES

## FEATURES

- All hardwood plywood with ultra-durable, UV-cured, satin finish base
- All base edging finished with tough Slate Gray PVC edging (no exposed plywood core)
- Ultra-rigid, steel uprights with durable finish
- Handrails and uprights are pre-assembled
- Each model is available in 30" and 36" widths
- Completely assembled with steel fasteners (no glued dowel pins)
- Reinforced bariatric capacity
- All steps are a full 10" deep
- Built for years of use



\* Available with optional Bus Step. See page 30.

# CLASSIC SERIES STAIRCASES

<b>4-5501-30</b>	Length	Width	Height
each step	10"	30"	6"
platform	24"	30"	24"
overall size	56"	34"	

### One-Sided Staircase

- Small corner staircase for limited space
- Requires a minimum 9' ceiling
- 500 lbs. load capacity

<b>4-5501-36</b>	Length	Width	Height
each step	10"	36"	6"
platform	24"	36"	24"
overall size	56"	40"	

### One-Sided Staircase

- Same as above with 36" width

Straight Position

<b>4-5500-30</b>	Length	Width	Height
1 step	10"	30"	6"
2 steps	10"	30"	4"
platform	24"	30"	12"
overall size	54"	34"	

### Small Up/Down Staircase

- Small straight staircase for limited space
- Patient can ascend and descend without turning
- Requires a minimum 8' ceiling
- 500 lbs. load capacity

<b>4-5500-36</b>	Length	Width	Height
1 step	10"	36"	6"
2 steps	10"	36"	4"
platform	24"	36"	12"
overall size	54"	40"	

### Small Up/Down Staircase

- Same as above with 36" width

Corner Position

<b>4-5080-30</b>	Length	Width	Height
3 steps	10"	30"	6"
5 steps	10"	30"	4"
platform	30"	30"	24"

### Convertible Staircase

- Converts from straight to corner arrangement
- 500 lbs. load capacity
- Requires a minimum 9' ceiling

<b>4-5080-36</b>	Length	Width	Height
3 steps	10"	36"	6"
5 steps	10"	36"	4"
platform	36"	36"	24"

### Convertible Staircase

- Same as above with 36" width