



SELECT SERIES STAIRCASES

FEATURES

- 1 1/2" diameter stainless steel handrails with end bumpers
- Fast easy height adjustments with self indexing spring action plungers and numbered height indicator strips (no holes to line up or count)
- Second knob removes any side-to-side play for a rigid upright
- Hardwood plywood base with durable, UV satin, topcoat finish
- Assembled with steel fasteners (No glued wood dowels)
- All edges are banded with tough, Slate Gray, PVC edge material
- 2 full width safety treads on each step
- Available in 30" and 36" widths
- 450 lbs. weight limit (both models)
- Reinforced bariatric load capacities



Adjust uprights by number strip on the uprights



Spring-loaded plungers make adjustments easy

4-6501-30	Length	Width	Height
platform	24"	30"	24"
each step	10"	30"	6"
overall size	54"	34"	

Adjustable One-Sided Staircase

- Adjusts in 1 1/2" increments to accommodate all ages, from children to adults
- Requires a minimum 9' ceiling

4-6501-36	Length	Width	Height
platform	24"	36"	24"
each step	10"	36"	6"
overall size	54"	40"	

Adjustable One-Sided Staircase

- Same as above with 36" width

4-6500-30	Length	Width	Height
platform	24"	30"	12"
2 steps	10"	30"	6"
3 steps	10"	30"	4"
overall size	54"	34"	

Adjustable Small Staircase

- Adjusts to all ages from children to adult in 1 1/2" increments
- Patient can ascend and descend without turning
- Requires a minimum 8' ceiling

4-6500-36	Length	Width	Height
platform	24"	36"	12"
2 steps	10"	36"	6"
3 steps	10"	36"	4"
overall size	54"	40"	

Adjustable Small Staircase

- Same as above with 36" width