

VALUE SERIES STAIRCASES

FEATURES

- Comfort grip moulded hardwood handrails
- Hardwood, plywood base with durable, satin, topcoat finish
- All stair edges are banded with tough Slate Gray, PVC edge material
- 2 full-width safety treads on each step
- Rigid grip handrail
- Fully enclosed sides helps give patient a secure feeling
- Assembled with steel fasteners (No glued dowel pins)
- 350 lbs. weight limit

	Length	Width	Height
4-7500	54"	34"	36"
platform	24"	30"	12"
1 step	10"	30"	6"
2 steps	10"	30"	4"

Small Straight Staircase with Full Sides

- Small straight staircase for limited space
- Patient can ascend and descend without turning
- Requires a minimum 8' ceiling

	Length	Width	Height
4-7501	56"	34"	58"
platform	24"	30"	24"
each step	10"	30"	6"

One-Sided Staircase with Full Sides

- Small corner staircase with three enclosed sides
- Requires a minimum 9' ceiling

101

Bus Step Option

- Creates a 12" high step when in place
- Available option for any staircase with a 6" rise step
- MUST BE ORDERED WITH STAIRCASE