

WORK CONDITIONING

	Width	Height	Depth
5-3078	30'	78"	20"

Heavy Duty Lift and Load Weight Shelf

- All hardwood, plywood construction
- Double thick uprights add durability
- All exposed edges are finished with bumper tee-moulding
- 3 adjustable shelves with easy lock-into-position stops
- 1" thick shelves are covered with durable, black, ABS plastic
- The complete unit is assembled with steel fasteners and special triple bolt steel corner brackets
- Upper shelf adjusts from 55" to 59"
- Middle shelf adjusts from 38" to 45"
- Lower shelf adjusts from 22" to 29"
- Maximum weight capacity per shelf 150 lbs.
- Must be attached to the wall for safety

ABS plastic on shelves

Double thick sides

	Width	Height	Depth
5-1515	15'	15"	15"

Square Combo Weight Box

- All hardwood, plywood construction
- Chrome center handle
- Approx. 16 lbs. empty, 100 lbs. capacity
- Angled top grip slots make it easy to carry

Triple bolted corners

← Triple Bolted Corners on each box

Bottom runners →

	Width	Height	Depth
5-1510	15'	15"	15"

Square Packing Weight Box with Lid

- All hardwood, plywood construction
- Approx. 21 lbs. empty, 100 lbs. capacity

	Width	Height	Depth
5-1915	19'	15"	15"

Rectangular Combo Weight Box

- All hardwood, plywood construction
- Chrome center handle
- Approx. 18 lbs. empty, 100 lbs. capacity
- Angled top grip slots make it easy to carry

	Width	Height	Depth
5-1717	17'	15"	17"

Large Storage Weight Box

- All hardwood, plywood construction
- Approx. 18 lbs. empty, 100 lbs. capacity
- Angled top grip slots make it easy to carry

	Length	Width	Height
5-2442	24'	48"	32"–45"

Heavy Duty Weight Sled

- 2 stainless steel push handles that are adjustable in height and depth
- All welded steel uprights with easy, spring loaded, locating indexing plungers
- Height is easily indicated by easy-to-read number indicator strip
- 4 removable chrome weight poles
- All hardwood, plywood weight box with metal triple bolted, 7 1/2" high, corner brackets
- Tapered self lubricating polyethylene sled runner
- Weight limit 150 lbs.
- Approximately 70 lbs. empty



Handrail depth adjustment

Indexing plunger and height indicator strip

5-0016

Disc Weight Set

- Black plate weights:
(4) 10 lbs.; (8) 5 lbs.;
(4) 2.5 lbs.

Flat

Vertical

	Width	Height	Depth
5-1924	25"	7 1/2"	16"

Multi Position Bolt Box

- Can be used flat or vertical
- Durable, ABS black, plastic nut and bolt shield
- 2 screwdrivers, 2 wrenches, 2 Allen wrenches plus 60 various nuts, bolts and washers are included

	Width	Height	Depth
5-3672	36'	62"	72"

Shovel Station

- 3 gravity flow shelves
- Folding side keeps sand in place
- 3 safety treads
- Tee-molded bumper edges
- All hardwood, plywood construction
- Plastic shovel and storage hook included

	Length	Width	Height
5-2323	23'	23"	8"

Extra Large Step Stool

- All plywood construction
- Safety tread top
- Bolt-through design