

Chemotherapy Spill Kit

*Protect and clean up spills involving cytotoxic therapy drugs**

This kit features Green-Z[®] solidifier to solidify and contain hazardous chemical spills for quick and easy clean-up and disposal with minimal exposure.

Item#	Description	Case Qty.	Case Wt.
49725	Chemotherapy Spill Kit (Poly Bag)	12	8 lbs.

*This kit may help you meet hazardous drug safety standard USP <800>. For more information on USP <800> compliance, visit readyfor800.com

Each Kit Contains:

- Nitrile Chemotherapy Gloves** (2 pairs)
- Chemotherapy Gown*** (1)
- Shoe Covers, non-slip (1 pair)
- Full-Face Safety Shield (1)
- 21g Green-Z[®] Fluid Solidifier (1)
- Scoop & Scraper (1)
- Yellow Chemotherapy Waste Bag (1)
- Wiper Pads (2)
- Twist Tie (1)
- Identification Tag (1)
- Instruction Card (1)
- Caution Sign (1)



Compliance Information

**Nitrile gloves tested for use with Chemotherapy drugs using ASTM D6978. (They are powder free, latex-free, DEHP free and silicone free)

***The Chemotherapy gown is compatible with OSHA recommendations, and ASHP and ONS guidelines for wearing a gown while preparing or administering chemotherapy drugs. The gown is disposable, lint-free, made of low permeability fabric, constructed with a solid front and back closure and is long sleeved with tight-fitting cuffs.